

Resona I9 Series

Innovation, in every facet Diagnostic Ultrasound System

Datasheet

Release 02.01.00

1 System Overview

Powered by the most revolutionary ZONE Sonography® Technology, I9's new ZST+ platform brings the ultrasound image quality to a higher level by zone acquisition and channel data processing.

It is intended for use in obstetrics, gynecology, abdominal, pediatric, small organ, musculo-skeletal, cardiac, vascular, urology, nerve, cephalic, thoracic/pleural, Laparoscopic, and Intra-operative exams.

1.1 Imaging Modes

- B-Mode
- THI and PSH (Phase Shift Harmonic Imaging)
- M-Mode/Color M-mode
- Free Xros M (Anatomical M-mode)
- Free Xros CM (Curved Anatomical M-mode)
- Color Doppler Imaging
- Power Doppler Imaging/Directional PDI
- PW (Pulsed Wave) Doppler
- CW (Continuous Wave) Doppler
- TDI (Tissue Doppler Imaging)
- Contrast Imaging
- Smart 3D (Freehand 3D)
- Real-time 4D
- iScape View (Panoramic Imaging)
- STE Imaging (Sound Touch Elastography)
- STQ Imaging (Sound Touch Quantification)
- Strain Elastography

1.2 Features

- B-Mode
- THI and PSH (Phase Shift Harmonic Imaging)
- M-Mode/Color M-mode
- Color Doppler Imaging
- Power Doppler Imaging and Directional PDI
- PW (Pulsed Wave) Doppler
- CW (Continuous Wave) Doppler
- Free Xros M
- Free Xros CM
- Glazing Flow
- iBeam (Spatial Compound Imaging)

- iClear (Speckle Suppression Imaging)
- iClear⁺
- iTouch (Auto Image Optimization)
- Echo Boost
- Zoom/iZoom (Full Screen Zoom)
- FCI (Frequency Compound Imaging)
- B steer
- ExFOV (Extended Field of View)
- HD Scope
- SSC (Sound Speed Compensation)
- Channel data processing
- Smart 3D
- Color 3D
- Real-time 4D
- STIC (Spatial-Temporal Image Correlation)
- iPage⁺ (Multi-Slice Imaging)
- SCV⁺ (Slice Contrast View)
- iLive
- Niche
- 3D-Print Format
- Smart Planes CNS
- Smart FLC
- Smart Planes FH
- Smart ICV
- Smart Face
- Smart-V (Smart Volume)
- Smart V Trace
- Smart Scene 3D
- Smart ERA
- IOTA
- Clinical Measurement Package
- Smart OB (Auto OB measurement)
- Smart NT (Auto NT measurement)
- Smart Fetal HR (Fetal Heart Rate)
- Smart HRI
- Smart Bladder
- Smart Hip
- Smart Trace
- Smart Calc
- CPP (Color Pixel Percentage)
- Smart Track
- Smart VTI
- Smart IVC
- Smart B-line
- HR Flow (High Resolution Flow)
- UMA (Ultra-Micro Angiography)
- IMT

- RIMT (RF-Data based IMT)
- R-VQS (RF-Data based Quantitative Analysis on Vessel Stiffness)
- Smart Pelvic
- Smart Breast
- Smart Thyroid
- IVF
- iScape View
- iNeedle (Needle Visualization Enhancement)
- V Flow
- Contrast Imaging
- Contrast Imaging QA (Quantitative Analysis)
- Volume CEUS
- LVO (Left Ventricular Opacification)
- Low MI Contrast
- CEUS Chrono-Parametric Mode
- TCMR
- STE Imaging (Sound Touch Elastography)
- STQ Imaging (Sound Touch Quantification)
- Strain Elastography
- High frame rate STE
- Endocavity STE
- USAT
- HRI⁺
- LTI
- Ultrasound Fusion Imaging
- Endocavity Fusion Imaging
- Fusion RESP
- ECG function
- AutoEF
- TDI (Include TVI, TVD, TVM, TEI)
- TDI QA (TDI Quantitative Analysis)
- TT QA (Tissue Tracking Quantitative Analysis)
- FH Tissue Tracking QA
- Stress Echo
- V-Mapping
- iScanHelper
- iWorks (Auto Workflow Protocol)
- DICOM
- MedSight
- MedTouch
- UltraAssist (Off-line software)
- UltraView (Off-line analysis software)

- Q-Path
- iStorage
- Touch gestures
- iVocal

2 Physical Specification

2.1 Dimension and Weight

The control panel and the monitor are in the lowest position.

- Configured with dual-wing floating arm and 23.8-inch monitor
 - Depth: 1020±20 mm
 - Width: 550±10 mm
 - Height: 1000±20 mm
- Weight: approx. less than 105kg±4kg (net weight, standard configuration but not including the transducer)

2.2 Electrical power

- Voltage: 100-240V~
- Frequency: 50/60 Hz
- Power consumption: Max. 650 VA

2.3 Operating Environment

- Ambient temperature: 0-40°C
- Relative humidity: 20%-85% (no condensation)
- Atmospheric pressure: 700hPa-1060hPa

2.4 Storage & Transportation Environment

- Ambient temperature: -20-55 °C
- Relative humidity: 20%-95% (no condensation)
- Atmospheric pressure: 700hPa-1060hPa

2.5 System Noise

26dB @25°C

3 User Interface

3.1 Monitor

- 23.8-inch high resolution color LED monitor
- Resolution: 1920x1080
- Viewing angle: 176 degrees
- Digital on screen display of brightness and contrast controls

- Automatic LED brightness
 - Tilt/Rotate independent adjustment
 - Tilt angle range: 20±5° degrees (backward), 85±5° degrees (forward)
 - Rotate angle range: 90±5° degrees (to left), 150±5° degrees (to right)
- 3.2 Multi-directional articulating monitor arm
- From left to right: 300±20mm
 - From front to back: arm: 300±20mm
 - From bottom to top: 150±20mm
- 3.3 Touch screen
- 15.6-inch high sensitivity anti-glare color touch screen
 - Resolution: 1920*1080
 - Digital brightness and contrast adjustment through preset
 - Viewing angle: 170 degrees
 - Angle adjustable range: 40 degrees
 - Support touch screen gestures
 - Support either hand writing or with gloves on
 - Movable 3D/4D tabs
 - Editable touch-screen buttons: long press to add, delete or move the buttons.
 - Digital TGC
 - Short-cut switch of latest used transducer & exams
- 3.4 Touch gestures
- Swipe down/up: display/remove projected image on touch screen
 - Swipe horizontally: page up/down or review images/cine loops one by one
 - Swipe from left edge to right: display hidden menu on projected image.
 - Image parameter adjustment.
 - Measurement on projected image on touch screen
 - Zoom in/out the projected image on touch screen
 - Rotate or erase on projected 3D/4D image on touch screen
 - 8 user defined gestures using two fingers for more functions, such as freeze, save, print, activate specific imaging modes, measurements, and some other special functions.
- 3.5 Floating control panel
- Brightness adjustable for the backlight of the whole control panel
 - Full-sized, backlit QWERTY keyboard
 - iConsole: intelligent control panel for clinical-exam specific layout and adaptive adjustment, 6 programmable E-ink keys for dynamic display of user-defined functions
 - Full-space floating control panel adjustment and can be fixed at any position (when centered in the trackball):
 - Left/right rotation: 80 degree for both left and right
 - Down/up adjustment: 300±20mm
 - Front/back adjustment: 200±20mm
- 3.6 Transducer port and holder
- Transducer ports with dust prevention: 5 active ports and 1 pencil transducer port
 - Support active transducer with indicator on
 - Transducer holder: 5, plus 1 dedicated endocavity transducer holder and 1 dedicated pencil transducer holder
- 3.7 Wheels
- Diameter: 125mm
 - When the central brake is configured: 3 castors for total lock, and 1 castor for direction lock and break.
 - When the central brake is not configured: 3 castors for total lock and break, and 1 castor for direction lock and break.
- 3.8 System boot-up
- Boot-up from shut-down: ≤ 30 sec
 - Boot-up from stand-by: < 5 sec
 - Shut-down: < 30 sec
- 3.9 Language support
- Support multi-language user interface, keyboard input, and user manuals.
 - User interface: English, Chinese, German, Spanish, French, Italian, Portuguese, Russian, Czech, Polish, Turkish, Finnish, Danish, Icelandic,

- Norwegian, Swedish, Hungarian, Serbian, Dutch, Lithuanian, Greek, Thai
- Keyboard input: English, Chinese, German, Spanish, French, Italian, Portuguese, Russian, Czech, Polish, Icelandic, Norwegian, Swedish, Finnish, Turkish, Danish, Hungarian, Serbian
- User manual: English, Chinese, German, Spanish, French, Italian, Portuguese, Russian, Polish, Turkish, Serbian, Norwegian, Danish, Swedish, Finnish, Dutch

* Not all items are listed in this part; For more detailed information, please refer to User Manual

3.10 Other Features

- DVD R/W driver
- Probe Container (For LAP13-4Cs)
- Anti-virus software
 - McAfee
 - Windows Defender
- Audio speakers: Stereo audio speakers
- Built-in Battery
 - Replaceable and rechargeable lithium battery.
 - Full battery lasts for no less than 24H in standby mode
 - Battery capacity indicators without powering on the system
 - Battery fully-recharged time: less than 6h
 - Continuous scanning time: more than 120 mins for 2H battery or 240 mins for 4H battery
- Built-in DVR
 - Built-in digital video recorder.
 - Screen display and voice information are recorded and stored in the built-in hard disk
 - Save space and is a useful tool for education and memory
 - Max storage length each time: 60 min

4 Technical Specifications

4.1 Innovative ZST⁺ platform

- Premium and innovative ultrasound platform, evolving with powerful

processing architecture and enhanced channel data processing based on ZONE Sonography®

- Windows Operating System
- Powerful Processing Architecture
- Advanced Acoustic Acquisition
- Total Recall Imaging (TRI)
- Enhanced Channel Data Processing
- Dynamic Pixel Focusing (DPF), digital variable aperture and dynamic apodization, A/D≥12 bit
- Up to 8,257,536 channels
- Parallel processing of multiple signals
- Line density per frame ≥ 512 ultrasound lines
- Sound Speed Compensation (SSC)
- ZONE Sonography® Technology

4.2 Transducer specification

- Single crystal and compound crystal etc.
- Supported transducer type
 - Curved array
 - Linear array
 - Phased array
 - Endocavity
 - 4D Volume
 - Bi-Plane
 - pencil transducer
 - Laparoscopic

5 Imaging Parameters

5.1 Advanced imaging technologies

- iBeam
Spatial Compound Imaging, permits use of multi-angle scanning to form a single image, so as to improve image contrast and resolution. Supports Spatial Compound Imaging with 3-level adjustment and up to 9-beam steer
- iClear
Speckle Suppression Imaging, Available in 2D, Contrast, and 3D/4D mode. 7-level adjustment.
- iTouch (Auto Image Optimization)
Automatic image parameters optimization of B, Color, PW, Contrast

- B-mode: gain, TGC, dehaze
 - Color: gain, color box position
 - Power: gain
 - PW: gain, scale, PRF, WF, SV size, SV position, steering angle
 - Contrast imaging: gain
- TSI

Multiple imaging conditions are available according to different tissue characteristics (general/muscle/fluid/fat)
- Smart Track
 - Available on linear transducers
 - Enable the function under Color/Power mode, the angle and the position of the ROI are adjusted automatically.
 - Enable the function under Color/Power mode. The angle and position of the ROI are automatically adjusted in real time, so as to save time from manual adjustment and improve the scanning efficiency.
- HD scope

By processing channel data multiply and retrospectively, HD Scope can improve the detail information and image contrast on specific area maximally.
- Echo Boost
 - Available in cardiac exam mode when using a phased array transducer
 - Improve the homogeneity through the whole field of view
 - Better improve the contrast display of the tissue
 - Better noise control in cardiac chambers and muscles
- iNeedle

It is used in needle-guided biopsy, which dynamically enhance the needle display in ultrasound image. Support dual-screen display of images before and after enhancement in real-time, and adaptive angle correction.

 - Needle direction: left, right
- B/iNeedle: on/off
- In-plane biopsy and Out-plane biopsy
- Zoom: Image magnification, supports spot zoom (write zoom) and pan zoom (read zoom). Magnification $\geq 10x$
- iZoom: Full-screen magnification by one click in real time, supports ≥ 2 magnification modes.
- Glazing Flow
 - An advanced image processing technology to demonstrate the 2D blood flow with 4D visualization for easy definition and clarity.
 - Available in Color/Power mode.
 - To stereoscopically display the blood flow.
- QSave
 - Quick save image parameter settings after image adjustment done
 - Support Save, Create, Restore
 - IP (Image Process): Quickly switch all the image parameters of the same image application by one-click
 - 3D/4D Preset Manager

The scenarios and subpresets can be renamed, restored, deleted, added, set to active, or moved

Show scenario and subpreset parameters

Provide multiple groups of preset 3D/4D parameters based on different application scenarios to quickly obtain expected image effect
- Auto Merge

Dual B image merge for linear array transducer
- Adjustable Acoustic output power
 - Display TIB, TIC, TIS in real time.
 - Can be adjusted in B, Color/Power, M, PW, CW, TDI modes etc.

5.2

B-mode

- Display formats: Single, Dual, Quad
- iClear/ iClear⁺: Off, 7 steps
- iBeam: Off, 3 steps
- iTouch: -12-12

- Dual Live: image compared by dual display in real time
 - Image quality
 - Supports fundamental frequency conversion ≥ 3 segments
 - Fundamental frequency: Pen/Gen/Res
 - Supports THI (Tissue Harmonic Imaging) and PSH (Phase Shift Harmonic Imaging), harmonic frequency conversion ≥ 3 segments
 - Harmonic frequency: HPen/HGen/HRes/HGen-FFR/HRes-FFR (dependent on transducer)
 - B steer: 5 levels, available on linear transducers
 - ExFOV
 - Available on transducers: convex, linear, endocavity, volume, phased, bi-planar, Laparoscopic.
 - Range: Off, 1-2 (dependent on transducer)
 - For liner transducer: after ExFOV is enabled, the image is displayed as a trapezoid. The maximum steer angle is 12°
 - For convex transducer: after ExFOV is enabled, the scanning angle will be extended.
 - Depth: 30 levels, 1.5-40cm (dependent on transducer)
 - Maximum display Depth: ≥ 40 cm
 - Maximum detection Depth of the convex transducer: ≥ 30 cm
 - Transmitting beam focusing: transmitting ≥ 8 segments
 - Supports FCI (Frequency Compound Imaging)
 - Frame rate (max): 1488 f/s
 - Acoustic output power: dependent on transducer
 - TGC: 8 segments on touch screen
 - LGC: 8 segments on touch screen
 - Dynamic range (visible and adjustable): 30-350 (dependent on probe/preset)
 - Gain: 0-100, 1/step
 - The Gain of B/M/D mode is visible and independently adjustable, ≥ 100 dB, 1dB/step
 - FOV Size: continuously adjustable
 - FOV Position: continuously adjustable (CB10-4s is not supported.)
 - Line density: L, M, H, UH
 - Persistence: 0-7 levels
 - Horizontal Scale: on/off
 - L/R flip and U/D flip: on/off
 - Rotation: 0° , 90° , 180° , 270°
 - TSI: general/muscle/fluid/fat
 - Gray Map: 8 types, grayscale: 256
 - Tint map: off, 8 types
 - Smooth: 0-6 levels
 - HD Scope: off, 1-3 levels
 - SSC (Sound Speed Compensation): on/off
 - Free view: adjustable (depends on transducers)
 - Dehaze: 0-30 levels
 - Ref Lines: on/off (under GYN and Pelvic Floor exam mode when using endocavity transducer)
 - V 1:1: on/off (available with linear transducer and under dual-split mode)
 - XL View: on/off
 - Edge Enhancement: 0-6
 - Echo Boost: on/off
 - Auto Merge: on/off
 - ZoneVue
- 5.3 THI and PSH
- Patent PSH technology, obtains purer harmonic, better contrast resolution, higher SNR, exceptional high frequency harmonic
 - iClear available
 - Image quality: HPen/HGen/HRes or HPen/HPen-FFR/HGen/HRes/HRes-FFR (depends on transducers)
 - Echo Boost: on/off
- 5.4 M-mode
- Display formats: V2:3, V3:2, V3:1, H2:3, FULL (V: vertical, H: horizontal)
 - Supports B/M mode dual display and M mode full screen display
 - Color M-mode available
 - Acoustic output power: same as B
 - Depth: same as B
 - Dynamic range: 30-180, 5/step

- Gain: 0-100, 1/step
- M sweep speed: 6 steps
- M soften: 0-4, 1/ step
- Tint map: off, 8 types
- Gray Map: 8 types
- Edge enhancement: 0-3, 1/ step

5.5 Color Doppler Imaging

- Dual live
- HR Flow: High Resolution Flow provides better image quality and flow sensitivity
- Image quality: Pen/Gen/Res (color), 1 level (HR Flow)
- Max velocity: 148.5 cm/s
- Steer: available on linear transducers
- Max frame rate: 279 f/s
- Acoustic output power: same as B mode
- Gain: 0-100, 2/step
- ROI size/position: adjustable
- Scale: max. 30 steps
- Baseline: -8 – 8, 1/step
- Wall filter: 8 steps, 5-433 Hz
- PRF: 0.1-15.4 kHz
- Packet size: 0-3, 1/ step
- Flow state: L/M/H
- Smooth: 0-6, 1/ step
- B/C align: on/off
- Priority: 0%-100%, 1%/step
- Color map: V0-V10; VV0-VV9
- Invert: on/off
- Auto Invert: on/off
- Persistence: 0-6, 1/ step
- Velocity tag: on/off
- Line density: L/M/H/UH
- iTouch: on/off
- Smart track: on/off
- ART Flow: Enhance the blood sensitivity and penetrability in time period
- Glazing flow: on/off, L/M/H

5.6 Power Doppler Imaging

- Dual live
- HR Flow: High Resolution Flow provides better image quality and sensitivity
- Support directional power Doppler

- Image quality: Pen/Gen/Res (Power), 1 level (HR Flow)
- Acoustic output power: same as B
- Dynamic range: 10-70, 5/step
- Gain: 0-100, 2/step
- ROI size/position: adjustable
- Steer: available on linear transducers)
- Scale: max. 30 steps
- Wall filter: 8 steps
- PRF: 0.1-15.4 kHz
- Packet size: 0-3, 1/ step
- Flow state: L/M/H
- Smooth: 0-6, 1/ step
- B/C align: on/off
- Priority: 0%-100%, 1%/step
- Color map: P0-P3, dP0-dP3
- Persistence: 0-6, 1/step
- Line density: L/M/H/UH
- Invert: on/off
- iTouch: on/off
- Smart track: on/off
- ART Flow: Enhance the blood sensitivity and penetration in certain time period
- Glazing flow: on/off, L/M/H

5.7 UMA (Ultra-Micro Angiography)

- The breakthrough of traditional Doppler imaging bottleneck. It realizes ultra-high spatial resolution and flow sensitivity for super subtle and slow flow imaging.
- Available under B, Color, Power, Contrast, and 3D mode.
- Supports submodes of cUMA, pUMA, sUMA.
- Supports a set of parameters for quick adjustment predefined for different organs in different scenes

Note: Other parameters are the same as those of the Color/Power modes

5.8 PW/CW Mode

- Display formats: V2:3, V3:2, V3:1, H2:3, FULL, Duplex/Triplex (PW only) (V: vertical, H: horizontal)
- Image quality: Pen/Gen/Res
- PW velocity:
 - max. 868.1 cm/s
 - min. 0.01 cm/s

- CW velocity:
 - max. 3900 cm/s
 - min. 0.01 cm/s
- Sample volume size: 0.5-30mm (PW only)
- Sample gate depth: continuously adjustable
- Baseline: 9 steps
- PW Steer: available on linear transducers
- Volume: 0%-100%, 2%/step
- PW PRF: 1.0-23.1 kHz
- CW PRF: 0.1-104.0 kHz
- Gain: 0-100, 2/step
- Dynamic range: 24-72, 2/step
- Sweep speed: 6 steps
- Wall filter
 - PW: 14-1000 Hz
 - CW: 7-1200 Hz
- Invert: on/off
- Auto invert: on/off
- Angle correction: -89~89 degrees, 1/step
- Quick angle: 0, -60, 60 degrees
- Gray map: 10 types
- Tint map: Off; 8 types
- Time/frequency resolution: 0-6, 1/step
- HPRF: On/Off
- Auto calc: on/off
- Auto calc cycle: 1, 2, 3, 4, 5
- Auto Calc Loop: on/off
- Trace Sensitivity: -3~3, 1/step
- Trace Smooth: -2, -1, off, 1, 2
- Trace area: above, below, all

5.9 Free Xros M

- Display formats: V2:3, V3:2, V 3:1, H2:3 (V: vertical, H: horizontal)
- Color Free Xros M available
- Up to 3 lines
- Display all lines
- Sweep speed: 6 steps
- M Tint map: off, 8 types
- Gray Map: 8 types

5.10 Free Xros CM

- Only available in TDI mode
- Display formats: V2:3, V3:2, V 3:1, H2:3 (V: vertical, H: horizontal)

- Sweep speed: 6 steps
- Tint map: off; 8 types
- Gray Map: 8 types
- Edit, undo, delete function for curved line

5.11 TDI Imaging

- Imaging modes: TVI, TEI, TVD, TVM
- Spectral Doppler frequency: ≥ 5
- Max frame rate: 1587 f/s
- Tissue Velocity/Energy Imaging (included in TDI option)
 - Available on transducers SC6-1s/SC8-2s/SP5-1s/P10-4s/P8-2s/P7-3Ts/P8-2Ts/P8-3Ts
 - Dual live: side by side displays B and B+TVI
 - PRF: 0.4-14.9 kHz
 - Acoustic output power: same as B mode
 - Gain: 0-100, 2/step
 - Dynamic range: 10-70, 5/step (TEI only)
 - ROI size/position: adjustable
 - Scale: max. 30 steps
 - Baseline: -8 – 8, 1/step (TVI only)
 - Wall filter: 8 steps
 - Packet size: 0-3, 1/step
 - Tissue state: L/M/H
 - Smooth: 0-6, 1/step
 - B/C wide: on/off
 - Priority: 0%-100%, 1%/step
 - Color map: 10 types (TVI), 8 types (TEI)
 - Invert: on/off
 - Persistence: 0-6, 1/step
 - Velocity tag: on/off (TVI only)
 - Line density: L/M/H/UH
 - Image quality: 2 levels
- Tissue Velocity Doppler (included in TDI option)
 - Display formats: V2:3, V3:2, V3:1, H2:3, FULL, Duplex/Triplex (V: vertical, H: horizontal)
 - Sample volume size: same as PW
 - Sample volume depth: continuously adjustable
 - Scale: 30 levels
 - PRF: 0.7-23.1 kHz

- Baseline: 9 levels
- Gain: 0-100, 2/step
- Dynamic range: 24-72, 2/step
- Sweep speed: 6 steps
- Wall filter: 10 steps
- Invert: on/off
- Angle correction: -89-89 degrees, 1/step
- Quick angle: 0, -60, 60 degrees
- Gray map: 10 types
- Tint map: Off; 8 types
- Image quality: 2 levels
- Time/frequency resolution: 0-6, 1/step
- Tissue Velocity Motion (included in TDI option)
 - Display formats: V2:3, V3:2, V 3:1, H2:3, FULL (V: vertical, H: horizontal)
 - Acoustic output power: same as B
 - Gain: 0-100, 2/step
 - M sweep speed: 6 steps
 - Smooth: 0-6, 1/step
 - Color Map: 10 types
 - Image quality: 2 levels
 - Persistence: 0-6, 1/step
 - Packet size: 0-3, 1/step
 - Priority: 0%-100%, 1%/step
 - Velocity tag: on/off
 - Tissue state: L/M/H

5.12 3D/4D

- Supports Freehand 3D (linear, convex, phased transducers)
- Supports real-time 4D and Color 3D (abdominal volume and endocavity volume transducers)
- General Imaging modes: Surface, Min, Max, X-ray
- Supports multi virtual light sources: Point, Parallel, Torch etc. and free combination of light sources is supported.
- Supports off-line processing of 3D/4D data, and adjusts imaging before storage of stored data.

5.13 Smart 3D

- Smart 3D Acquisition preparation:

- 3D/4D Scenario setting: Routine
- Acquisition Methods: Rocked, Linear
- Reset VOI: On/Off
- Flip VOI: On/Off
- Angle: 10-80°
- Distance: 10-200mm
- Acquiring Time: 1.0s-20.0s

VR:

- 3D/4D Scenario: Routine-Surf., iLive Gen., iLive Transp., Skeleton
- Reset: All, Orientation, Curve
- VOI: On/Off/Fixed
- Active quadrant: A, B, C, VR
- VR orientation: 0°, 90°, 180°, 270°
- Flip: flip VR
- Sync: synchronize VR with selected plane
- Orientation Assist: On/Off
- Threshold: 0-100%, 1%/step
- Opacity: 0-100%, 5%/step
- Brightness: 0-100%, 2%/step
- Contrast: 0-100%, 2%/step
- Smooth: 0-10, 1/step
- Depth VR: Off/Black/Cyan/Blue/Rose
- Tint: off; 8 types
- Degree: 10-80°
- Distance: 10-200mm

MPR:

- Active quadrant: A, B, C
- Gray Map: 1-8
- Brightness: 0-100%, 2%/step
- Contrast: 0-100%, 2%/step
- iClear: Off; 7 types
- Tint: Off; 8 types
- Thickness: 0-30mm

Adv.:

- Direction: Up/Down, Left/Right, Front/Back, Down/Up, Right/Left, Back/Front
- VR Refine: Off; 7 steps
- Surface enhance: 0-7, 1/step
- MagiClean: Off/Low/Mid/High/Max
- Inversion: On/Off
- A3:1: On/Off
- Move light: On/Off

- Degree: 10-80°
- Distance: 10-200mm
- Main render: Surface, Max, Min, X Ray, iLive
- Sub render: Surface, Max, Min, X Ray
- Mix: Set the mix ratio of the two render modes
- Tool
 - Edit:
 - Rubber: On/Off
 - Eraser Diameter: 8-80, 1/step
 - Cut (area selection): Polygon, Contour, Rectangle
 - Undo: Undo, Undo All
 - 3D Layout:
 - Niche Views: Inner, Outer
 - 3Slice
 - Active Quadrant: A, B, C, Niche/3Slice
 - Auto rotation:
 - Position: Set Start/Set End
 - Direction: Left/Right, Up/Down
 - Step: 1-15°
 - Quick Angle: 30-180°
 - Rotation control: play, single loop, loop
 - Save AVI to USB

5.14 4D

- Available on all volume transducers
- Static 3D and real time 4D
- Acquisition preparation:
 - 4D frame rate: max. 80 vps
 - 3D/4D Scenario setting: Smart Scene3D (Spine, Brain, Long Bone, Face, Endometrium, Pelvic), Routine, iLive Pro, Bone, Tissue (not all scenarios are listed)
 - Refresh: On/Off
 - Angle
 - Quality: low1, low2, mid, high1, high2
- VR:
 - 3D/4D Scenario: Smart Scene 3D, Routine, iLive Pro, Bone, Tissue, Routine (not all scenarios are listed)
 - 3D iClear: Off; 7 steps

- Face⁺: Off, 3 steps
- Auto Play: Stop, x1, x2, x3, x1/2, x1/3
- Frame: Select a frame
- Other parameters are the same as Smart 3D

MPR

- 3D iClear: Off; 7 types
- Other parameters are the same as Smart 3D

Adv.:

- Threshold: 0%-100%
- Opacity: 0%-100%
- Brightness: 0%-100%
- Contrast: 0%-100%
- 3D iClear: Off, 1-7
- Smooth: 0-10
- Depth VR: OFF, Black, Cyan, Blue, Rose
- Tint: Off, 1-8
- Shading: 0-10
- Hyaline: 0%-100%
- Grad View
- Other parameters are the same as Smart 3D

Tool

The parameters are the same as Smart 3D

- 3D Reference Point

- Enable the operator to define one or more reference points on MPRs, which are then projected to VR image; helpful for the operator to better understand the corresponding spatial relations of VR image and MPRs
- Display: Point, H line, V line
- Delete All
- Hide All

- 3D Print

- Quality: Low, Mid, High
- Generate Mesh
- File Format: .stl, .obj, .ply, .3mf, .off
- Save mesh to USB
- Threshold: 0%-100%
- 3D iClear: Off, 1-7
- Smooth: 0-10

5.15 Color 3D

Provides more stereoscopic blood flow signals. It is mainly used in blood perfusion tissues in some complex space to make blood flow observation visually. 3D imaging in color and power modes is supported.

5.16 Niche

Compiles the 3 MPRs together according to their relative positions, to provide a much clearer interior anatomical structure for diagnosis

5.17 Smart Scene 3D

An innovative technique of automatic scenario-oriented volume scan for extremely easy, efficient and accurate exams. It is capable of automatically identifying tissue characteristics, such as fetal brain, fetal face, fetal bone, fetal spine

5.18 STIC

Spatio-Temporal Imaging Correlation, based on the movements of the fetus, rebuilt and show the anatomical structure within a physical movement by using the interconnection between the time and the space, for better diagnosis in 4D function.

5.19 iPage⁺

- Multiple tomographic parallel slices imaging, to display volume image with slices, for better display the spatial relationship of tissue and lesion.
- Displays ≥ 25 images at different depths on the screen at the same time, and the slice spacing is adjustable (0.5 mm-10 mm).

5.20 SCV⁺

Slice Contrast View Plus, includes SCV imaging and CMPR functions. SCV imaging can reduce speckle noise and improve contrast resolution as well as enhance signal-noise ratio, which helps in discovering diffuse pathology in organs. CMPR is to obtain a curved sectional image by dissecting the 4D image, so as to observe the curved tissue structure. In clinical application, it is often used to observe the curved anatomy in stretching such as fetal spine.

5.21 iLive

- An advanced rendering mode for realistic volume imaging display
- Brings a better imaging experience by adding lighting rendering effect to the traditional way, allowing human tissue texture to be revealed more clearly.
- Supports Hyaline function.

5.22 Smart Planes CNS

- A professional screening software of the fetal central nervous system. Automatically detect the standard CNS scanning planes, and calculate the anatomical parameters for each planes. It is used to assist doctors in diagnosis of fetal craniocerebral diseases in a more efficient and standardized way.
- Automatically acquire multiple standard cranial sections and acquire ≥ 4 commonly used measurement indicators
- Automatically display the standard planes: TCP, TTP, MSP and TVP
- Auto comment supported: A(anterior), P(posterior), L(Left), R(Right), U(Up), D(down), CSP, T, CH, CV, CM, LV on TCP, TTP, MSP and TVP
- Auto measurement supported:
 - TCD and Cist Maga (CM) on section TCP;
 - BPD, OFD and HC on section TTP;
 - LVW on section TVP
- Support editing measurement results
- Hide/show measurement results
- Support comment and bodymark on sectional plane

5.23 Smart ICV

- Smart intracranial volume, an advanced tool to detect fetal cranial tissue, provides automatic calculation of fetal intra-cranial volume for advanced fetal CNS study.
- Supports modifying the contour by trace, adaptive trace, control point editing, and recalculate the volume.

5.24 Smart Planes FH

- Detect automatically left ventricular outflow tract view, right ventricular outflow tract view, LAV-DA view, 3VV-T view and stomach bubble view
- Automatically acquire ≥ 6 standard fetal heart sections.

5.25 Smart Face

Allows to recognize fetal face and remove the shading obstacle data automatically, then display the face in an optimal viewing angle. At the same time, you can adjust the display direction of the fetal face by one click, and support forward/reverse rubber erasing.

5.26 Smart FLC (Smart Follicle)

- Automatically measures and calculate the number and size of follicles in the image area, and evaluate follicles according to the follicle size.
- Automatically segments the anechoic structure is by one click, and displays the anechoic structures at different positions and sizes in different colors.
- Automatically measures the follicle diameter, length at the X-axis, Y-axis, Z-axis, as well as the average value and volume of the three axes.

5.27 Smart-V (Smart Volume)

Fast volume calculation tools to calculate the volume of tissue structure or lesions

- Smart-V ROI: Manual ROI on A, B, C plane separately
- Smart-V Trace: Manual trace on A, B, C plane separately
- Smart-V Vocal/Smart-V Parallel: Trace contours on each slice separately

5.28 Smart ERA

- A fully automated endometrial receptivity analysis tool, Enables endometrium receptivity assessment with automated workflow.
- Supports automatic measurement of the segmented results.

5.29 RIMT (RF-Data based IMT)

- Available in single/dual B carotid exam mode

- Side: left/right
- Calculation of 6 RIMT values, RIMT average value, SD and ROI W
- Report operation:
- Data deleting
- RIMT trend graphic viewing
- Preview

5.30 iScape View

- Acquisition method: B and Power
- Supports speed indicator
- Actual size: on/off
- Fit size: on/off
- Ruler: on/off
- Tint map: off; 8 types
- Rotation: 0~355 degrees, 5/step

5.31 V Flow (Vector Flow)

- A novel approach for vascular hemodynamic analysis, using color coded vector arrows to display flow velocity magnitude and direction,
- Available on L9-3s/L14-3Ws/L13-3Ns in real-time B/Color carotid exam mode
- Quality: 1-8, 1/step
- V Flow map: 4 levels
- Gain: 0-8, 1/step
- Persistence: 0-5, 1/step
- Edge smooth: Low/Mid/High
- Wall filter: 8 steps
- Brightness: 0-100%, 5%/step
- Colorbar scale: 50-100cm/s, 10cm/s/step
- Arrow life time: 10-50, 1/step
- Arrow size: S/M/L/XL/XXL
- Arrow density: 1%-10%, 1%/step
- Packet Size: 0-3, 1/step
- Acquire Time: 0.5, 0.6, 0.8, 1.0, 1.2, 1.5s

5.32 Contrast Imaging*

- Contrast imaging technology, which provides exceptional contrast agent detecting capability, not only extracts second harmonic, but also non-linear fundamental signals
- Available on C6-2Gs/C11-3s/SC6-1s/SC8-2s/C4-1s/C9-3Ts/CB10-4s/V11-3Hs/V11-3HBs/L9-3s/L14-3Ws/L20-5s/L16-4Hs/7LT4s/L12-3RCs/L13-3Ns/SP5-1s/LAP13-4Cs/DE11-3Ws/SD8-1s/ELC13-4s

- Micro Flow Enhancement (MFE) available
- Timer1: on/off
- Timer2: on/off
- Pro capture: captures prospective image less than 480s preset table
- Retro capture: captures retrospective image less than 120s preset table
- Dual live: side by side displays tissue image and contrast image
- MFE period: 0.1s, 0.2s, 0.4s, 0.6s, 0.8s, 1.0s, MAX
- Destruct: instantly destroy contrast bubbles
- Destruct AP: -30~0 dB, 0.3/step
- Destruct time: 500-2000 ms,75/step
- iClear: off; 7 steps
- Mix: mix contrast image with tissue image
- Mix map: 7 types, available when Mix mode is active
- Persistence: 8 steps
- Dynamic range: same as B mode
- Gray map: 8 types
- Tint map: off; 8 types
- Supports U/D Flip and L/R Flip
- Rotation: same as B mode
- CEUS Position: on/off
- Line density: L/M/H/UH
- FOV: on/off
- FOV size/position: continuously adjustable
- ExFov: off, 1-2, 1/step
- Gain: 0-100, 1/step
- iTouch: on/off, -8~8, 2/step
- Image quality: 3 levels
- Smooth: 0-6, 1/step
- Enhance: on/off
- Markline: on/off
- LGC: 8 points
- ZoneVue
- Tissue Gain: 0-100,1/step
- *This ultrasound series is designed for compatibility with commercially available ultrasound contrast agents. Because the availability of these agents is subject to government regulation and approval, product features intended for use with these agents may

not be commercially marketed nor made available before the contrast agent is cleared for use. Contrast related product features are enabled only on systems for delivery to an authorized country or region of use. Mindray medical systems makes no claims concerning the safety or effectiveness of contrast agents.

5.33 Contrast Imaging QA

- Support Time-Intensity Curve analysis
- Table display: display data in table
- Up to 8 ROIs
- Delete all
- Delete current
- Fit curve
- Raw curve
- Motion tracking: Reduce the effect of tissue movement
- X scale: 1-5, 1/step
- Export: export current data as CSV format file

5.34 LVO

- Only available in cardiac exam mode
- Dedicated left ventricle contrast imaging tool

5.35 Low MI Contrast

- Only available in cardiac exam mode
- Enhances echo reflection by using contrast agent to perform myocardial analysis based on echocardiography technique

5.36 Volume CEUS

- Available on DE11-3Ws/SD8-1s
- Timer1: on/off
- Timer2: on/off
- Capture 3D image

5.37 CEUS Chrono-Parametric Mode

- To intuitively visualize the difference between the target lesion and surrounding tissue, CCPM depicts temporal information as a color coded map superimposed on the CEUS image.
- Available under 2D and 4D mode.
- Background: set the background image for the color coding of arriving time.

- Motion Tracking: enable or disable motion tracking compensation.
- Transparency: adjust the transparency for the color coding of arriving time.
- Threshold: adjust the threshold of the contrast agent arriving time.
- Color Map: used during color coding of the contrast agent arriving time in the contrast image.

5.38 TCMR

Tissue-Contrast Mix Rendering, a mixed rendering mode with both contrast perfusion and tissue information for intuitive display of the spatial relationship and helpful for easy positioning and lesion diagnosis.

5.39 STE Imaging (Sound Touch Elastography Imaging)

- Available on C6-2Gs/SC6-1s/SC8-2s/C4-1s/L9-3s/L14-3Ws/L20-5s/L13-3Ns
- Endocavity STE
 - CB10-4s, V11-3Hs, DE11-3Ws, ELC13-4s
 - 2D real-time STE of the endocavity transducer. Mainly used for urology and gynecology examinations.
- Display Format: V1:1, H1:1, FULL
- Invert: on/off
- HQ Elasto: on/off
- E Quality: Pen, Gen, Res
- Elas.Metric: E, Cs, G
- Scale: 30 levels
- Opacity: 0-5, 1/step
- Map: 3 types
- ROI Width/Height: continuously random adjustable
- ROI Center Depth: continuously adjustable
- Depth: same as B mode
- iLayering: on/off
- Filter: 0, 1
- RLB View: on/off
- M-STB Index: on/off
- M-STB Sensibility: 0~4, 1/step
- Smooth: 0~2, 1/step
- Persistence: 0~2, 1/step
- RLB Map: on/off, RLB, RLB&E, RLB&B&E

- Map Position: 0%~100%, 5%/step
- E bar: Mean, Max, Min, SD
- E Avg: off, 8 levels
- Select/Bad: on/off
- Lesion: off, 1~10
- FS Mode: on/off

5.40 High frame rate STE

To increase the frame rate of STE function.

5.41 Strain Elastography

- Available on V11-3Hs/V11-3HBs/L9-3s/L14-3Ws/L20-5s/L16-4Hs/7LT4s/L12-3RCs/L13-3Ns/DE11-3Ws/ELC13-4s
- Support strain, strain ratio and strain histogram measurement
- Unique shell analysis function
- Stress compensation technology reduces deeper tissue artifacts, obtain more uniform stress throughout whole field
- Stress indicator: supports frame by frame stress indication.
- Map: 6 types
- Smooth: 0-5, 1/step
- Opacity: 0-5, 1/step
- ROI Width/height: continuously adjustable
- Invert: on/off
- Display Format: V1:1, H1:1, FULL
- Strain mode: 0~1, 1/step
- Dynamic Range: 0~5, 1/step
- Map Position: 0%~100%, 5%/step

5.42 STQ Imaging (Sound Touch Quantification Imaging)

- Available on C6-2Gs/SC6-1s/SC8-2s/C4-1s/CB10-4s/V11-3Hs/L9-3s/L14-3Ws/L20-5s/L13-3Ns/DE11-3Ws/ELC13-4s
- ROI Adjustment: adjust the ROI fixed size
- Elas.Metric: E, Cs, G
- The square height of the elasto curve represents the average value of the elasto metric for current frame.
- E bar: Mean, Max, Min, SD
- M-STB Index: On/Off
- M-STB Sensibility: 0-4, 1/step
- Filter: 0, 1
- Smooth: 0-2, 1/step

- Persistence: 0-2, 1/step
- Map Position: 0~100%, 5%/step
- Lesion: off, 1~10
- Scale: 0-9, 1/step
- E Avg: off, 8 levels
- HQElasto: on/off
- FS Mode: on/off

5.43 Fatty Liver Lab

Only the SC6-1s transducer supports this function

- USAT
UltraSound ATtenuation analysis enables quantitative fatty liver assessment by measuring the attenuation coefficient.
- HRI⁺
HepatoRenal Index Plus enables quantitative fatty liver assessment by measuring the echo intensity ratio between the liver parenchyma and the renal cortex based on the RF data.
- LTI
Liver Texture Index enables quantitative fatty liver assessment by statistically analyzing the size and density of scatters of the liver parenchyma
- Parameters of USAT, HRI⁺, and LTI
Note: the parameter items may vary based on the specific fatty liver lab mode selected
 - Scan Time: Single, 1s, 2s, 3s, MAX
 - H-ROI Size: 10mmx10mm~26mmX26mm(2mm/step), 30mmx30mm
 - M-STB Index: On/Off
 - Filter: 0, 1, 2
 - M-STB Sensi.: 0~4,1/step
 - H1:1: On/Off
 - Full: On/Off
 - E bar: Mean, Max, Min, SD (the value may varies based on the specific fatty liver lab mode selected)
 - RLB Map: On/Off
 - RLB View: On/Off
 - Opacity: 0~9, 1/step
 - Map: F1~F4, 1/step

- Att.Metric: dB/cm/MHz, dB/m/MHz, dB/m
- Scale: 0~10, 1/step
- Frame Average: 1, 3, 5, 7, 10 (the value may varies based on the specific fatty liver lab mode selected)
- Select/Bad: 1~10, 1/step

5.44 Ultrasound Fusion Imaging

- Available on C6-2Gs/SC6-1s/SC8-2s/C4-1s/L14-3Ws/L13-3Ns/SP5-1s in B/ Color/ Power/ contrast imaging mode (non-cardiac)
- Single window display
- Fusing CT/MR/PET/freehand volume data with the ultrasonic image
- CT/MR/PET data reconstruction for 3D displaying
- Tracking system: on/off
- Fusion ratio: -1~1, 0.1/step
- Axis rotation: 0° ~ 360° based on X-axis, Y-axis or z-axis in increment of 2°.
- ROI Offset X: -630~630
- ROI Offset Y: -566~566
- Window W/L: 1/step. Adjust the CT/MR/PET/freehand brightness and the contrast by changing the width and the level.
- Reset Window W/L
- Reset CT/MR
- Display marks
- Respiration curve: on/off
- Respiration Range: 0/1/2
- View Type: Axial, Coronal, Sagittal
- Quick Translation
- Registration
- Mark on Volume Data/Mark on Fusion Image
- Support general measurement
- Support adding comment and bodymark

5.45 Endocavity Fusion Imaging

Mainly used for endocavity transducer, fusing real-time ultrasound and CT/MR images.

5.46 Fusion RESP

By automatically detecting the patient's respiratory motion signal, provide the respiratory motion compensation technology, reduce the fusion distortion caused by the patient's respiration, and provide more accurate fusion imaging.

5.47 AutoEF

- Automatic Ejection Fraction Measurement, it provides an efficient way to detect the left ventricle and calculate the EF.
- Output EDV/ESV/EF/SV/CO by Simpson method
- Activated with or without ECG
- Adjustment for the border of endocardium by single point or multi points
- Adjust Frame
- Layout: Dual/ Single
- Diastole FR
- Systole FR
- Volume curve: on/off

5.48 TDI QA

- Dedicated quantification tool for TDI velocity, strain, strain rate analysis
- Ellipse ROI, Standard ROI
- Up to 8 of ROI
- ROI tracking: tracking ROI along with cardiac movement
- Delete all
- Delete current
- Smooth: 1-7, 1/step
- X scale: 1-5, 1/step
- Std.Height: 1.5-50 mm
- Std.Width: 1.5-50 mm
- Std.Angle: -89-90 degrees
- Export: export current data as CSV format file

5.49 TT QA

- Available on SP5-1s/P10-4s/P8-2s in adult cardiac, cardiac-difficult (car-penetration), and pediatric cardiac exam modes.
- Tissue tracking quantitative analysis
- Mandatory ECG connection before TT QA cine acquisition

- Six views for analysis: ALAX, A4C, A2C, PSAXB, PSAXM, PSAXAP
- Reload: reload cine again for new study
- Edit: modify trace points
- Start tracking
- Accept & compute: start tracking myocardium movement when user accept trace result
- Display effect: 0/1; at 0, tracking in dots; at 1, tracking in velocity vector arrow
- Trace method: 3 point or manual for ALAX, A4C, A2C; manual for PSAXB, PSAXM, PSAXAP
- Bull's eye: trace result in bull's eye model
- LGC: available
- Valve's open and close time index: MVC, MVC; AVC, AVO, MVO
- Data export: export data in CSV file
- Cycle: ECG triggered cardiac cycle recognition for analysis; adjustable
- Auto play: stop, X1/10, X1/5, X1/4, X1/3, X1/2, X1, X2, X3
- Thickness: 1-30mm, 1mm/step; adjust trace thickness
- Track point: 20-40, 1/step
- Parameter: Volume, Speed, Displacement, L Strain, L Strain R, T Strain, T Strain R, Area, R Strain, R Strain R, C Strain, C Strain R, C Rotation, C Rotation R
- Smooth: 0-4, 1/step
- Trace method: 3 point, manual
- Tracking cycles: 1, 3
- Blue's Eye view
- Curve Display
- Torsion & Torsion Rate Curve
- LGC adjustment

5.50 Fetal Heart TT QA

- Fetal Heart Tissue Tracking with Quantitative Analysis, Real-time track the continuous motion of fetal myocardium, display the mechanical changes of each segment in the longitudinal, circumferential and radial directions, offering accurate and effective evaluation of myocardial movement of fetal heart

- Available Under Fetal Echo related exam modes
- 5.51 Stress Echo
- Available on SP5-1s/P10-4s/P8-2s in cardiac mode
 - 14 factory protocols
 - User-defined protocols
 - ECG triggered acquisition, display, selection, comparison, evaluation and archiving of multiple cardiac loops during various stages of a stress echo examination
 - Customized stages: up to 7 views per stage, and up to 12 stages per study
 - View: standard views (PLAX, SAB, PSAX, SAA, A4C, A2C, ALAX), and customized views
 - Image acquisition
 - R-wave trigger
 - Acquire mode: Manual ROI or full screen
 - Ability to acquire frames or clips in B-mode, M-mode, Color, PW and TDI
 - Image selection
Attach the images with view annotation label (PLAX, SAB, PSAX, SAA, A4C, A2C, ALAX, and customized views)
 - Review
Automatically adjust to the number of images user defined
 - Wall Motion Scoring
 - ASE 16 (with score 4-7), or ASE 17(with score 4-7)
 - Graphical display of scoring (Normal, Hyperkinetic, Severely Hyperkinetic, Akinetic, Dyskinetic)
 - LV volume measurement
Measurement of LV Volume in all phases of cardiac cycle
 - Report
Reporting for both Wall Motion Scoring and LV volume measurement
- 5.52 Smart IVC
- Automatic Inferior Vena Cava trace and calculation, automatically trace the IVC diameter change, and calculate the CI, DI and IVC Variation, helping for volume status assessment and guide the fluid therapy.
- 5.53 Smart VTI
- Automatic Velocity Time Integral measurement, automatically trace the PW spectrum of LVOT, and obtain VTI, SV, CO and SVV, for rapid assessment of volume responsiveness.
- 5.54 Smart B-line
- Used to detect B lines of the lung in B mode
 - Acquisition method: single B in real-time or freeze mode
 - Scanning areas: 6 Zones, 8 Zones, 12 Zones
 - Auto Calc
 - OverView
 - Image and diagnosis comments
- 5.55 R-VQS
- RF-data Based Quantitative Analysis on Vessel Stiffness
 - Track movements of the upper and lower vessel walls automatically
 - Displacement and Vessel diameter display in the result window.
 - Motion curve of vessel walls display under the image in real time.
- 5.56 Smart Pelvic
- Enter smart pelvic in 2D or 3D/4D scanning mode.
 - Set Rest and Valsalva frame
 - Measure automatically
- 5.57 Smart Fetal HR (Fetal Heart Rate)
- Measure the fetal heart rate automatically on B/M mode
- 5.58 Smart Breast
- Automated breast mass BI-RADs classification and reporting tool, which provides automated breast ultrasound lesion classification and reporting, and it can be used as breast ultrasound quality control solution with high-efficiency.
 - Standardized scanning procedure based on the lesions, and Standardized scanning procedure guide

- Automatic lesion identification, detection, measurement, comment, report generation, etc.
 - Multi-lesion& multi-plane management and analysis
- 5.59 Smart Thyroid
- Automated thyroid ultrasound analysis and reporting tool, for efficiency and standardization of thyroid nodule diagnosis.
 - Standardized scanning procedure based on the lesions, and Standardized scanning procedure guide
 - Automatic lesion identification, detection, measurement, comment, report generation, etc.
 - Multi-lesion& multi-plane management and analysis
- 5.60 IVF
- In Vitro Fertilization, includes dedicated IVF image parameter/ measurement/comment/bodymark / report
 - The uterus and follicle growth curve can be displayed in the IVF report.
 - Data of IVF history exams can be checked in the IVF report.
 - The user-defined analysis model is supported for reproductive function evaluation.
- 5.61 iScanHelper
- An ultrasonic tutorial and assistant package integrated in ultrasound devices.
 - Tutorial function as a guidance to show basic scanning skill with graphic of transducer position, schematic of anatomy and example clinical image
 - Supports Abdomen, gynecological, urological, obstetrical, Small Parts and nerve block area.
- 5.62 iCompare
- Allow to compare real-time ultrasound imaging to images from iStation; Helpful to easily evaluate and follow up the progression of disease, treatment effect monitoring
- 5.63 V-Mapping
- Manually edit the vascular body mark.
 - Manually edit the body mark through the touch screen to Intuitively display the lesion position.
- 6 Cine Review and Raw Data Processing
- 6.1 Cine review
- Available in 2B, Color, Contrast, 4D modes.
 - Supports editing and clipping during cine review. Frame by frame manual cine loop review or auto playback with variable speed
 - Save and edit the clipped and edited cine images for many times
 - Maximum cine memory up to 44107 frames (B storage server) or 134.05s (M storage server) (depends on the mode)
 - Maximum 4D cine memory: 10721 volumes (DE11-3Ws)
 - Retrospective storage (1-120s pre-settable) and prospective storage (1-480s pre-settable)
 - Frame compare: displays one cine in dual format and allows frame by frame compare side by side
 - Cine compare: compare cines which are saved in same imaging mode
 - Jump to first and jump to last: one keystroke go to first or last frame in the cine
- 6.2 Raw data processing
- B-mode
 - TGC
 - Gain
 - Dynamic range
 - Gray map
 - Tint map
 - iClear
 - L/R Flip
 - U/D Flip
 - Rotation
 - iTouch
 - LGC
 - Dual live

- Auto Merge
- H Scale
- Echo Boost
- B/iNeedle
- Smooth
- Zoom
- Ref Lines
- Dehaze
- V1:1
- XL View
- Edge Enhance
- M-mode
 - Gain
 - Speed
 - Dynamic Range
 - Gray Map
 - Tint Map
 - Display format
- Color
 - Gain
 - Baseline
 - Smooth
 - Color map
 - Priority
 - Dual Live
 - Invert
 - Velocity tag
 - Glazing flow
- PW
 - Gain
 - Baseline
 - Volume
 - Angle
 - Speed
 - Dynamic range
 - Gray map
 - Tint Map
 - Display format
 - Invert
 - WF
 - T/F Res

measurement packages, see the Appendix.

7.1 Automatic Measurement and Analysis

- AutoCalc
 - PS
 - ED
 - MD
 - PPG
 - TAMAX
 - Vol Flow(TAMAX)
 - TAMEAN
 - Vol Flow(TAMEAN)
 - Vas Diam
 - Vas Area
 - DT
 - MPG
 - MMPG
 - VTI
 - AT
 - S/D
 - D/S
 - PI
 - RI
 - PV
 - HR
- IMT (Intima-Media Thickness Measurement)

Automatic detection, measurement and analysis of IMT when ROI is set
- Smart OB
 - Auto measurement for OB, a special tool for easy OB scan, and greatly reduce time and increase productivity
 - Support BPD, HC, OFD, FL, AC, HUM. More than 6 automatic measurement items are supported.
- Smart NT

Automatic detection and measurement of fetal nuchal translucency, which makes NT evaluation more efficient and accurate.
- Smart HIP

Automatic hip measurements and Graf classification of neonatal and pediatric hip for easy screening of DDH (Developmental Dysplasia of the Hip).

7 Measurement/Analysis and Report*

NOTICE:

For general measurement, automatic measurement, and clinical

- **Smart HRI**
Hepato-Renal Index, which is provided automated liver steatosis assessment by automatically calculating brightness scale ratio between the liver with renal cortex in B mode.
- **Smart Trace**
A smart tool for easy and precise boundary trace
Automatically identify and trace lesions, as well as measure the length, area and circumference of the long axis and short axis within the defined ROI.
- **CPP**
Measure the blood flow signal distribution in the defined ROI under the Color/Power mode, and calculate the color pixel percentage within the defined ROI.
- **Smart Bladder**
Auto measure three diameters and calculate the bladder volume
- **Smart Calc**
Automatic trace, measurement and calculation tool

7.2 User-defined Measurement

Supports user-defined measurement calculations and studies

7.3 Report

- Specific report template by application
- Editable value in report
- Images selectable
- Able to Export as PDF file
- Presetable hospital information
- Anatomy information for vascular and OB report
- Editing though iReport
- User-defined report template
- Selecting report modules
- Adding/removing measurement items from the report
- Changing report layout
- Load/save comment
- Viewing history reports
- Preview and printing reports

7.4 Comments/Bodymark

- Supports text input and arrow

- Voice annotation: record voice as annotation for images and cine
- Support freehand marking on touch screen
- Adjustable text size and arrow size
- Supports home position
- Covers various application
- More than 250 bodymarks for versatile application
- User customizable Comments/ Bodymark

7.5 iWorks

- Auto workflow protocol, can be combined with a standardized exam protocol to facilitate a more efficient and standard ultrasound exam in the clinic with guidance.
 - Automatically add comments, body marks and switch the image mode according to the protocol.
 - Templates are user configurable
 - Template import and export are available
 - Functions: pause, stop, replace, repeat, skip, insert single step, return and continue, steps in thumbnail
 - iNSert another template during the iWorks process.
 - iWorks setup mode: B; B/B (Dual Live); Dual B/B; Color; B/Color (Dual Live); Power; B/Power (Dual Live); B + PW; Color + PW; Power + PW; B + CW; Color + CW; B+M; B+TVI; TVI+TVD; iScape View.
 - iWorks setup annotation: support up to 2 annotations, location and font size are configurable.
 - iWorks setup bodymark: select existing library, and transducer indicator is pre-settable
 - iWorks setup measurement: select existing measurement library
- * Not all measurements are listed in this part; For more detailed information, please refer to User Manual

8 Exam Storage and Management

8.1 Exam Management

- iStation workstation dedicated for patient information management
- Workstation dedicated for image management
- Workstation dedicated for report and report template management
- Workstation dedicated for exam management
 - Patient exam query/retrieve
 - Support review of current and past exam
 - New exam, Activate exam, Continue exam functions, End exam are available
 - Support measurements and calculations on archived exam and images

8.2 Exam Storage

- Supports local hard disk storage: 1TB HDD (Hard Disk Drive) and 128G SSD (Solid State Drive)
- Direct digital storage of single frame and cine files.
- Export images in BMP/JPG/TIFF/DCM/AVI/MP4/WMV format
- Supports back-end storage, export and backup of image data, real-time exam, without affecting exam operation.
- Supports backup/send to USB devices
- Supports backup/send to DVD-RW media
- Supports DICOM 3.0
- Supports network storage (iStorage). Based on the TCP/IP protocol, the ultrasound image and report can be directly transmitted to the PC device.
- Supports print service.

9 Connectivity

9.1 Ethernet Network Connection

- Cable connection
- Wireless connection

9.2 DICOM 3.0

- DICOM Basic

- Verify (SCU, SCP)
- Print
- Store
- Storage Commitment
- Media Exchange
- DICOM Worklist
- DICOM Query/Retrieve
- DICOM Modality Performed Procedure Step - MPPS
- DICOM OB/GYN structure report
- DICOM Cardiac structure report
- DICOM Vascular structure report
- DICOM Breast structure report
- DICOM Abdomen structure report
- DICOM Small Parts structure report
- DICOM Pediatric structure report
- DICOM Urology structure report

9.3 UltraView

Off-line analysis workstation, PC-based ultrasound image analysis software. Supports post-processing and more measurement analysis of ultrasound image off-line.

- Components:
 - DICOM Basic
 - DICOM Query/Retrieve
 - DICOM OB/GYN SR
 - DICOM Cardiac SR
 - DICOM Vascular SR
 - DICOM Breast SR
 - TDI QA
 - Contrast Imaging QA
 - Tissue Tracking QA
 - iPage
 - Niche
 - SCV (Slice Contrast View)
 - iLive
 - Smart OB
 - Smart NT
 - IVF
 - AutoEF
 - IMT
 - Stress Echo
 - Ultrasound Fusion Imaging
 - Smart Planes CNS

9.4 UltraAssist

Direct network storage tool between ultrasound system and personal computer

- The ultrasound system supports sending reports or patient information to the PC via iStorage.
- The PC supports import custom report templates and test items to the ultrasound system via iReport and iMeasurement.

9.5 MedSight

- An interactive App that lets you transfer clinical images straight from Mindray Ultrasound system to a smart device, such as mobile phone or tablet PC
- Transfer images or clips from system to mobile terminal through Wi-Fi
- Support both iOS and Android system
 - Android (4.0 and above)
 - iOS (7.0 and above)
 - For iOS powered smart device: DICOM is mandatory
 - For Android powered smart device: DICOM not necessary

9.6 MedTouch

- Connect Ultrasound machine to smart devices based on Android and iOS system, such as tablet PC or mobile phone. Remote control of Ultrasound machine, review of patient information, and tutorial software iScanHelper study on smart devices
- Supports remote control of the Ultrasound machine, to adjust parameter (Gain, Depth, Freeze/Unfreeze, Image Storage, Zoom, iTouch, mode switching among B, Color, Dual Live).
- Support Android and iOS powered smart devices
 - Android 4.0 and above
 - iOS 7.0 and above
 - DICOM not necessary

10 Transducers

10.1 Curved Array

- C6-2Gs
 - Application: Obstetrics, Gynecology, Abdominal, Urology
 - Advanced Function: iScape View, Free Xros M, Color M, Smart 3D, Contrast Imaging (Gynecology, Abdominal), STE (Abdominal), STQ (Abdominal), Ultrasound Fusion Imaging
 - Bandwidth: 1.2-6.0 MHz
 - Number of Elements: 128
 - FOV (max): 94°
 - Extended FOV: 106°
 - Convex Radius: 20 mm
 - Depth: 4.0-40 cm
 - Physical Footprint: 37.6 mm × 19 mm
 - Footprint: 31.5 mm × 11.2 mm
 - B-mode Frequencies: 1.2-3.8, 1.7-5.2, 2.0-6.0 MHz
 - Harmonic Frequencies: 3.2, 4.0, 5.0 MHz
 - Color Frequencies: 2.0, 2.5, 3.0, 3.5 (HR Flow) MHz
 - PW Frequencies: 2.0, 2.5, 3.0 MHz
 - Biopsy Guide: NGB-024, multi angle, reusable; LPUBKG81 (disposable)
- C11-3s
 - Application: Abdominal, Vascular, Cardiac, Small Organ, Pediatric, Cephalic
 - Advanced Function: iScape View, Free Xros M, Color M, Smart 3D, Contrast Imaging (Abdominal, Vascular)
 - Bandwidth: 2.6-12.8 MHz
 - Number of Elements: 128
 - FOV (max): 101°
 - Extended FOV: 113°
 - Convex Radius: 15 mm
 - Depth: 1.5-35 cm
 - Physical Footprint: 32.8 mm × 25 mm
 - Footprint: 27.4 mm × 8.4 mm

- B-mode Frequencies: 2.6-6.5, 3.2-7.9, 4.7-12.8 MHz
- Harmonic Frequencies: 6.0, 7.0, 8.0 MHz
- Color Frequencies: 4.4, 5.0, 5.7, 5.0 (HR Flow) MHz
- PW Frequencies: 4.4, 5.0, 5.7 MHz
- Biopsy Guide: NGB-018, multi angle, reusable
- SC6-1s
 - Application: Obstetrics, Gynecology, Abdominal, Small Organ, Musculo-skeletal, Vascular, Urology, Nerve, Thoracic/pleural
 - Advanced Function: iScape View, Free Xros M, Color M, Smart 3D, TDI (Obstetrics), Contrast Imaging (Gynecology, Abdominal, Vascular, Urology), STE (Abdominal), STQ (Abdominal), Ultrasound Fusion Imaging
 - Bandwidth: 1.2-6.0 MHz
 - Number of Elements: 192
 - FOV (max): 60°
 - Extended FOV: 72°
 - Convex Radius: 60 mm
 - Depth: 4.0-40.0 cm
 - Physical Footprint: 65.1 mm × 16.4 mm
 - Footprint: 64.9 mm × 16.2 mm
 - B-mode Frequencies: 1.2-3.8, 1.7-5.2, 2.0-6.0 MHz
 - Harmonic Frequencies: 4.0, 5.0, 6.0 MHz
 - Color Frequencies: 2.0, 2.5, 3.0, 3.3 (HR Flow) MHz
 - PW Frequencies: 2.0, 2.5, 3.0 MHz
 - Biopsy Guide: NGB-022, multi angle, reusable; LPUBKG60 (disposable); CIVCO 658-004-1/2/3/4/5 (disposable)
- SC8-2s
 - Application: Obstetrics, Gynecology, Abdominal, Urology, Vascular
 - Advanced Function: iScape View, Free Xros M, Color M, Smart 3D, TDI (Obstetrics), Contrast Imaging (Gynecology, Abdominal, Urology, Vascular), Ultrasound Fusion Imaging, STE (Abdominal), STQ (Abdominal)
 - Bandwidth: 1.8-8.2 MHz
 - Number of Element: 192
 - FOV (max): 76°
 - Extended FOV: 88°
 - Convex Radius: 40 mm
 - Depth: 4.0-40.0 cm
 - Physical Footprint: 26.3 mm × 66.9 mm
 - Footprint: 15 mm × 52 mm
 - B-mode Frequencies: 1.8-5.4, 2.3-6.8, 2.8-8.2 MHz
 - Harmonic Frequencies: 4.0, 5.5, 6.0 MHz
 - Color Frequencies: 3.0, 3.5, 4.0, 3.5 (HR Flow) MHz
 - PW Frequencies: 3.0, 3.5, 3.8 MHz
 - Biopsy Guide: NGB-029, multi angle, reusable
- C4-1s
 - Application: Gynecology, Obstetrics, Abdominal, Small Organ, Vascular, Urology, Thoracic/pleural
 - Advanced Function: iScape View, Free Xros M, Color M, Smart 3D, Contrast Imaging (Gynecology, Abdominal, Vascular, Urology), STE (Abdominal), STQ (Abdominal), Ultrasound Fusion Imaging
 - Bandwidth: 1.0-4.5 MHz
 - Number of Elements: 64
 - FOV (max): 56°
 - Extended FOV: 81°
 - Convex Radius: 30 mm
 - Depth: 4.0-40.0 cm
 - Physical Footprint: 42.4 mm × 25.4 mm
 - Footprint: 34.1 mm × 17.5 mm
 - B-mode Frequencies: 1.0-3.0, 1.6-3.3, 2.3-4.5 MHz
 - Harmonic Frequencies: 2.4, 3.0, 3.4 MHz
 - Color Frequencies: 1.6, 1.9, 2.3, 2.5 (HR Flow) MHz
 - PW Frequencies: 1.6, 1.9, 2.5 MHz
 - Biopsy Guide: NGB-036, multi angle, reusable; CIVCO 698-013

- (disposable); CIVCO 698-019 (disposable)
- **CB10-4s**
 - Application: Urology
 - Advanced Function: iScape View, Free Xros M, Color M, Smart 3D, Endocavity STE, STQ, Contrast Imaging
 - Bandwidth: 2.6-12.8 MHz
 - Number of Elements: 128
 - FOV (max): 165°
 - Extended FOV: 205°
 - Convex Radius: 9 mm
 - Depth: 1.5-28.0 cm
 - Physical Footprint: 22.5 mm × 20.1 mm
 - Footprint: 20.1 mm × 9 mm
 - B-mode Frequencies: 2.6-6.5, 3.2-7.9, 4.7-12.8 MHz
 - Harmonic Frequencies: 7.0, 8.0, 9.0 MHz
 - Color Frequencies: 4.4, 5.0, 5.0, 5.0 (HR Flow) MHz
 - PW Frequencies: 4.4, 5.0, 5.7 MHz
 - Biopsy Guide: NGB-004, single angle, reusable; CIVCO 610-543 (disposable); CIVCO 610-1274 (disposable)
- **C9-3Ts**
 - Application: Abdominal, Obstetrics, Musculo-skeletal, Vascular, Small Organ, Intra-operative, Pediatric
 - Advanced Function: iScape View, Free Xros M, Color M, Smart 3D, Contrast Imaging (Abdominal)
 - Bandwidth: 2.6-9.0 MHz
 - Number of Elements: 128
 - FOV (max): 68°
 - Extended FOV: 80°
 - Convex Radius: 33 mm
 - Depth: 4.0-40.0 cm
 - Physical Footprint: 46 mm × 37.5 mm
 - Footprint: 42 mm × 10.5 mm
 - B-mode Frequencies: 2.6-4.8, 3.6-6.4, 4.5-7.2, 5.5-9.0 MHz
 - Harmonic Frequencies: 4.5, 5.0, 6.0 MHz
- **V11-3Hs**
 - Color Frequencies: 3.0, 3.3, 3.5, 3.3 (HR Flow) MHz
 - PW Frequencies: 3.0, 3.3, 3.8 MHz
 - Biopsy Guide: not available
- **V11-3Hs**
 - Application: Obstetrics, Gynecology, Urology
 - Advanced Function: iScape View, Free Xros M, Color M, Smart 3D, Contrast Imaging (Gynecology, Urology), Strain Elastography (Gynecology, Urology), Endocavity STE (Gynecology, Urology), STQ (Gynecology, Urology), Endocavity Fusion Imaging (Urology)
 - Bandwidth: 3.0-11.0 MHz
 - Number of Elements: 192
 - FOV (max): 170°
 - Extended FOV: 210°
 - Convex Radius: 11 mm
 - Depth: 1.5-28.0 cm
 - Physical Footprint: 24.9 mm × 21.8 mm
 - Footprint: 24 mm × 9 mm
 - B-mode Frequencies: 3.0-7.0, 4.0-9.0, 5.0-11.0 MHz
 - Harmonic Frequencies: 8.0, 9.0, 10.0 MHz
 - Color Frequencies: 4.4, 5.0, 5.7, 5.5 (HR Flow) MHz
 - PW Frequencies: 4.5, 5.0, 5.5 MHz
 - Biopsy Guide: NGB-025, single angle, reusable; CIVCO 610-543 (disposable); CIVCO 610-1274 (disposable)
- **V11-3HBs**
 - Application: Obstetrics, Gynecology, Urology
 - Advanced Function: iScape View, Free Xros M, Color M, Smart 3D, Contrast Imaging (Gynecology, Urology), Strain Elastography (Gynecology, Urology)
 - Bandwidth: 3.0-11.0 MHz
 - Number of Elements: 192
 - FOV (max): 180°
 - Extended FOV: 240°
 - Convex Radius: 11 mm
 - Depth: 1.5-28.0 cm

- Physical Footprint: 23.4 mm × 19.2 mm
- Footprint: 23.8 mm × 9.1 mm
- B-mode Frequencies: 3.0-7.0, 4.0-9.0, 5.0-11.0 MHz
- Harmonic Frequencies: 8.0, 9.0, 10 MHz
- Color Frequencies: 4.4, 5.0, 5.7, 5.5 (HR Flow) MHz
- PW Frequencies: 4.5, 5.0, 5.5 MHz
- Biopsy Guide: NGB-048, single angle, reusable

10.2 Volume

- DE11-3Ws
 - Application: Obstetrics, Gynecology, Urology
 - Advanced Function: iScape View, Free Xros M, Color M, 3D/4D, Contrast Imaging (Gynecology, Urology), Volume CEUS (Gynecology, Urology), Strain Elastography (Gynecology, Urology), Endocavity STE (Gynecology, Urology), STQ (Gynecology, Urology)
 - Bandwidth: 2.0-9.0 MHz
 - Number of Elements: 192
 - FOV (max): 162°
 - Extended FOV: 187°
 - Volume Sweep Angle (max): 120°
 - Convex Radius: 10 mm
 - Depth: 4.0-40.0 cm
 - Physical Footprint: 24 mm × 24 mm
 - Footprint: 24 mm × 24 mm
 - B-mode Frequencies: 2.0-6.0, 2.8-8.2, 3.0-9.0 MHz
 - Harmonic Frequencies: 4.0, 5.0, 6.0 MHz
 - Color Frequencies: 4.4, 5.0, 5.7, 5.0 (HR Flow) MHz
 - PW Frequencies: 4.4, 5.0, 5.7 MHz
 - Biopsy Guide: NGB-047, single angle, reusable
- SD8-1s
 - Application: Obstetrics, Gynecology, Abdominal
 - Advanced Function: iScape View, Free Xros M, Color M, 3D/4D, Contrast Imaging (Gynecology,

Abdominal), Volume CEUS (Gynecology, Abdominal)

- Bandwidth: 1.8-8.2 MHz
- Number of Elements: 192
- FOV (max): 66°
- Extended FOV: 91°
- Volume Sweep Angle (max): 85°
- Convex Radius: 45 mm
- Depth: 4.0-40.0 cm
- Physical Footprint: 75.7 mm × 52.6 mm
- Footprint: 54.5 mm × 14.9 mm
- B-mode Frequencies: 1.8-5.4, 2..3-6.8, 2.8-8.2 MHz
- Harmonic Frequencies: 4.0, 5.0, 5.5, 6.0, 6.3 MHz
- Color Frequencies: 3.0, 3.5, 4.0, 3.5 (HR Flow) MHz
- PW Frequencies: 3.0, 3.5, 3.8 MHz
- Biopsy Guide: NGB-039, multi angle, reusable

10.3 Linear Array

- L9-3s
 - Application: Abdominal, Pediatric, Small Organ, Musculo-skeletal, Vascular, Nerve, Obstetrics
 - Advanced Function: CW, iScape View, Free Xros M, Smart 3D, V Flow (Vascular), Contrast Imaging (Abdominal, Pediatric, Small Organ, Musculo-skeletal, Vascular), Strain Elastography (Small Organ, Musculo-skeletal, Vascular), STE (Small Organ, Musculo-skeletal, Abdominal, Vascular), STQ (Small Organ, Musculo-skeletal, Abdominal, Vascular)
 - Bandwidth: 2.5-9.0 MHz
 - Number of Elements: 192
 - FOV (max): 4.37 cm
 - Extended FOV
 - Ohter: 20°
 - OB1, NT: 30°
 - Steered Angle
 - B: -12°, -6°, 0°, 6°, 12°
 - C/PW: -30°-30°
 - Depth: 1.5-35.0 cm
 - Physical Footprint: 62 mm × 22 mm

- Footprint: 48 mm × 11 mm
- B-mode Frequencies: 2.5-7.0, 3.4-8.2, 3.6-9.0 MHz
- Harmonic Frequencies: 5.0, 6.0, 7.0 MHz
- Color Frequencies: 3.0, 3.6, 5.0, 4.0 (HR Flow) MHz
- PW Frequencies: 3.0, 3.6, 5.0 MHz
- Biopsy Guide: NGB-034, multi angle, reusable
- L14-3Ws
 - Application: Small Organ, Musculo-skeletal, Vascular, Abdominal, Pediatric, Thoracic/ Pleural, Nerve
 - Advanced Function: iScape View, Free Xros M, Smart 3D, V Flow (Vascular), Strain Elastography (Small Organ, Musculo-skeletal, Vascular), STE (Small Organ, Musculo-skeletal, Vascular), STQ (Small Organ, Musculo-skeletal, Vascular), Contrast Imaging (Small Organ, Musculo-skeletal, Vascular, Abdominal, Pediatric), Ultrasound Fusion Imaging
 - Bandwidth: 3.0-14.0 MHz
 - Number of Elements: 256
 - FOV (max): 5.08 cm
 - Extended FOV: 20°
 - Steered Angle
B: -12°, -6°, 0°, 6°, 12°
C/PW: -30°-30°
 - Depth: 1.5-35.0 cm
 - Physical Footprint: 66.8 mm × 25.5 mm
 - Footprint: 55.5 mm × 8.2 mm
 - B-mode Frequencies: 3.0-9.0, 5.0-12.0, 6.0-14.0 MHz
 - Harmonic Frequencies: 10.0, 11.0, 12.0 MHz
 - Color Frequencies: 5.0, 6.2, 7.2, 8.3 (HR Flow) MHz
 - PW Frequencies: 5.0, 6.2, 7.2 MHz
 - Biopsy Guide: NGB-054, multi angle, reusable
- L20-5s
 - Application: Abdominal, Small Organ, Musculo-skeletal, Vascular, Nerve
- L16-4Hs
 - Application: Musculo-skeletal, Nerve, Small Organ, Vascular, Pediatric, Intra-operative
 - Advanced Function: iScape View, Free Xros M, Smart 3D, Strain Elastography (Small Organ, Musculo-skeletal, Vascular), Contrast Imaging (Small Organ, Musculo-skeletal, Vascular)
 - Bandwidth: 3.5-16.0 MHz
 - Number of Elements: 128
 - FOV (max): 2.53 cm
 - Extended FOV: 20°
 - Steered Angle
B: -12°, -6°, 0°, 6°, 12°
C/PW: -30°-30°
 - Depth: 1.5-35.0 cm
- L14-3Ws (continued)
 - Advanced Function: CW, iScape View, Free Xros M, Smart 3D, Strain Elastography (Small Organ, Musculo-skeletal), STE (Small Organ, Musculo-skeletal, Vascular), STQ (Small Organ, Musculo-skeletal, Vascular), Contrast Imaging (Small Organ)
 - Bandwidth: 6.0-23.0 MHz
 - Number of Elements: 192
 - FOV (max): 2.85cm
 - Extended FOV: 20°
 - Steered Angle
B: -12°, -6°, 0°, 6°, 12°
C/PW: -20°-20°
 - Depth: 1.5-29.0 cm
 - Physical Footprint: 42.23 mm × 22.10 mm
 - Footprint: 31.5 mm × 4.5 mm
 - B-mode Frequencies: 6.0-13.0, 9.0-16.6, 12.5-23.0 MHz
 - Harmonic Frequencies: 14.0, 16.0, 18.0 MHz
 - Color Frequencies: 9.0, 11.0, 13.0, 13.0 (HR Flow) MHz
 - PW Frequencies: 8.3, 10.0, 12.5 MHz
 - CW: 10.0 MHz
 - Biopsy Guide: not available

- Physical Footprint: 11.5 mm x 38 mm (slant width)/34.8 mm (straight width)
- Footprint: 28.7 mm × 5.5 mm
- B-mode Frequencies: 3.5-9.2, 7.6-12.8, 9.6-16.0 MHz
- Harmonic Frequencies: 10.0, 12.0, 14.0 MHz
- Color Frequencies: 6.2, 7.2, 8.3, 4.0 (HR Flow) MHz
- PW Frequencies: 5.0, 6.2, 8.3 MHz
- Biopsy Guide: not available
- 7LT4s
 - Application: Small Organ, Musculo-skeletal, Vascular, Abdominal, Pediatric, Intra-operative, Thoracic/Pleural
 - Advanced Function: iScape View, Free Xros M, Smart 3D, Strain Elastography (Small Organ, Musculo-skeletal, Vascular), Contrast Imaging (Abdominal)
 - Bandwidth: 3.5-13.5 MHz
 - Number of Elements 128
 - FOV (max): 3.98 cm
 - Extended FOV: 10°
 - Steered Angle
B: -12°, -6°, 0°, 6°, 12°
C: -15°-15°
PW: -30°-30°
 - Depth: 1.5-35.0 cm
 - Physical Footprint: 49.5 mm × 14.4 mm
 - Footprint: 45 mm × 9.0 mm
 - B-mode Frequencies: 3.5-8.3, 4.4-9.2, 5.6-13.5 MHz
 - Harmonic Frequencies: 8.0, 9.0, 10.0 MHz
 - Color Frequencies: 5.0, 6.5, 8.0, 8.0 (HR Flow) MHz
 - PW Frequencies: 5.0, 6.5, 8.0 MHz
 - Biopsy Guide: NGB-010, multi angle, reusable
- L12-3RCs
 - Application: Small Organ, Vascular, Musculo-skeletal, Abdominal, Pediatric, Thoracic/Pleural, Nerve
 - Advanced Function: iScape View, Free Xros M, Smart 3D, Contrast Imaging (Small Organ), Strain Elastography (Small Organ)
 - Bandwidth: 3.0-11.0 MHz
 - Number of Elements: 192
 - FOV (max): 3.80cm
 - Extended FOV: 20°
 - Steered Angle
B: -12°, -6°, 0°, 6°, 12°
C/PW: -30°-30°
 - Depth: 1.5-35.0 cm
 - Physical Footprint: 55.6 mm × 22 mm
 - Footprint: 43.5 mm × 8.2 mm
 - B-mode Frequencies: 3.0-8.3, 4.4-9.2, 5.6-11.0 MHz
 - Harmonic Frequencies: 7.0, 8.0, 9.0 MHz
 - Color Frequencies: 4.4, 5.0, 7.2, 6.0 (HR Flow) MHz
 - PW Frequencies: 4.2, 5.0, 7.2 MHz
 - Biopsy Guide: NGB-043 (within the plane), multi angle, reusable; NGB-044 (outside the plane), multi angle, reusable
- L13-3Ns
 - Application: Small Organ, Musculo-skeletal, Vascular, Abdominal, Pediatric, Thoracic/Pleural, Nerve
 - Advanced Function: iScape View, Free Xros M, Smart 3D, V Flow (Vascular), Strain Elastography (Small Organ, Musculo-skeletal, Vascular), STE (Small Organ, Musculo-skeletal, Vascular), STQ (Small Organ, Musculo-skeletal, Vascular), Contrast Imaging (Small Organ, Musculo-skeletal, Vascular, Abdominal, Pediatric), Ultrasound Fusion Imaging
 - Bandwidth: 3.0-13.0 MHz
 - Number of Elements: 192
 - FOV (max): 3.80 cm
 - Extended FOV: 20°
 - Steered Angle
B: -12°, -6°, 0°, 6°, 12°
C/PW: -30°-30°
 - Depth: 1.5-35.0 cm

- Physical Footprint: 56.8 mm × 21.2 mm
- Footprint: 43.5 mm × 8.2 mm
- B-mode Frequencies: 3.0-9.6, 5.4-11.5, 6.6-13.0 MHz
- Harmonic Frequencies: 9.0, 10.0, 10.5, 11.0 MHz
- Color Frequencies: 4.2, 5.0, 6.2, 7.2 (R Flow) MHz
- PW Frequencies: 4.2, 5.0, 6.2 MHz
- Biopsy Guide: NGB-053, multi angle, reusable

10.4 Phased Array

- SP5-1s

- Application: Abdominal, Cardiac, Vascular, Cephalic, Thoracic/ Pleural
- Advanced Function: CW, iScape View, Free Xros M, Free Xros CM, Color M, Smart 3D, TDI (Cardiac), LVO (Cardiac), Low MI Contrast (Cardiac), Stress Echo (Cardiac), Tissue Tracking QA (Cardiac), Ultrasound Fusion Imaging
- Bandwidth: 1.5-4.5 MHz
- Number of Elements: 80
- FOV (max): 90°
- Extended FOV: 90°
- Depth: 2.0-38.0 cm
- Physical Footprint: 38.2 mm × 30.5 mm
- Footprint: 23.4 mm × 15.2 mm
- B-mode Frequencies: 1.5-2.5, 2.5-3.5, 3.5-4.5 MHz
- Harmonic Frequencies: 3.0, 3.5, 4.0 MHz
- Color Frequencies: 2.0, 2.3, 2.5, 2.5 (HR Flow) MHz
TDI: 3.0, 3.8 MHz
- PW Frequencies: 2.0, 2.3, 2.5 MHz
TDI: 2.5, 4.0 MHz
- CW: 2.0 MHz
- Biopsy Guide: NGB-011, multi angle, reusable

- P10-4s

- Application: Abdominal, Cardiac, Pediatric, Nerve, Cephalic
- Advanced Function: CW, iScape View, Free Xros M, Free Xros CM,

Color M, Smart 3D, TDI (Cardiac), Stress Echo (Cardiac), Tissue Tracking QA (Cardiac)

- Bandwidth: 3.0-11.4 MHz
- Number of Elements: 128
- FOV (max): 90°
- Extended FOV: 90°
- Depth: 2.0-16.5 cm
- Physical Footprint: 15.1 mm × 10.2 mm
- Footprint: 15 mm × 9.1 mm
- B-mode Frequencies: 3.0-6.8, 3.8-10.2, 4.6-11.4 MHz
- Harmonic Frequencies: 7.5, 8.0, 9.0 MHz
- Color Frequencies: 4.0, 5.0, 5.7, 6.3 (HR Flow) MHz
TDI: 5.7, 6.2 MHz
- PW Frequencies: 4.0, 5.0, 5.7 MHz
TDI: 5.0, 5.7 MHz
- CW: 5.0 MHz
- Biopsy Guide: not available

- P8-2s

- Application: Abdominal, Cardiac, Nerve, Pediatric, Cephalic
- Advanced Function: CW, iScape View, Free Xros M, Free Xros CM, Color M, Smart 3D, TDI (Cardiac), Stress Echo (Cardiac), Tissue Tracking QA (Cardiac)
- Bandwidth: 2.3-8.0 MHz
- Number of Elements: 96
- FOV (max): 90°
- Extended FOV: 90°
- Depth: 2.0-38.0 cm
- Physical Footprint: 30.5 mm × 23.2 mm
- Footprint: 19.5 mm × 11 mm
- B-mode Frequencies: 2.3-5.4, 2.8-7.4, 4.2-8.0 MHz
- Harmonic Frequencies: 5.0, 6.0, 7.0 MHz
- Color Frequencies: 2.7, 3.3, 4.0, 4.0 (HR Flow) MHz
TDI: 5.0, 6.2 MHz
- PW Frequencies: 2.7, 3.3, 4.0 MHz
TDI: 5.0, 6.2 MHz
- CW: 2.5 MHz
- Biopsy Guide: not available

- P7-3Ts
 - Application: Cardiac
 - Advanced Function: CW, Free Xros M, Free Xros CM, Color M, TDI (Cardiac)
 - Bandwidth: 2.3-7.2 MHz
 - Number of Elements: 64
 - FOV (max): 90°
 - Extended FOV: 90°
 - Depth: 2.0-38.0 cm
 - Physical Footprint: 14 mm × 12 mm
 - Footprint: 12.2 mm × 12.2 mm
 - B-mode Frequencies: 2.3-5.4, 2.8-6.4, 3.3-7.2 MHz
 - Harmonic Frequencies: 5.0, 6.0, 7.0 MHz
 - Color Frequencies: 2.7, 3.3, 4.0, 4.0 (HR Flow) MHz
TDI: 5.0, 6.2 MHz
 - PW Frequencies: 2.7, 3.3, 4.0 MHz
 - Biopsy Guide: not available
 - P8-2Ts
 - Application: Cardiac
 - Advanced Function: CW, Free Xros M, Free Xros CM, Color M, TDI (Cardiac)
 - Bandwidth: 2.3-7.2 MHz
 - Number of Elements: 64
 - FOV (max): 90°
 - Extended FOV: 90°
 - Depth: 2.0-38.0 cm
 - Physical Footprint: 14 mm × 12 mm
 - Footprint: 10.6 mm × 10.6 mm
 - B-mode Frequencies: 2.3-5.4, 2.8-6.4, 3.3-7.2 MHz
 - Harmonic Frequencies: 5.0, 6.0, 7.0 MHz
 - Color Frequencies: 2.7, 3.3, 4.0, 4.0 (HR Flow) MHz
TDI: 5.0, 6.2 MHz
 - PW Frequencies: 2.7, 3.3, 4.0 MHz
TDI: 2.7, 5.0 MHz
 - CW: 2.5 MHz
 - Biopsy Guide: not available
 - P8-3Ts
 - Application: Cardiac
 - Advanced Function: CW, Free Xros M, Free Xros CM, Color M, TDI (Cardiac)
 - Bandwidth: 2.3-7.2 MHz
 - Number of Elements: 48
 - FOV (max): 90°
 - Extended FOV: 90°
 - Depth: 2.0-38.0 cm
 - Physical Footprint: 10.7 mm × 7.9 mm
 - B-mode Frequencies: 2.3-5.4, 2.8-6.4, 3.3-7.2 MHz
 - Harmonic Frequencies: 6.0, 6.5, 7.0 MHz
 - Color Frequencies: 3.3, 3.8, 4.4, 4.0 (HR Flow) MHz
TDI: 5.0, 6.2 MHz
 - PW Frequencies: 3.3, 3.8, 4.4 MHz
TDI: 5.0, 6.2 MHz
 - CW: 2.5 MHz
 - Biopsy Guide: not available
- 10.5 Bi-Plane
- ELC13-4s (Convex & Linear)
- Application: Urology, Gynecology
 - Advanced Function: iScape View, Free Xros M, Color M, Smart 3D, Strain Elastography, Endocavity STE, STQ, Contrast Imaging, Endocavity Fusion Imaging
 - Bandwidth
 - Convex: 3.5-9.5 MHz
 - Linear: 4.8-11.0 MHz
 - Number of Elements: 192
 - FOV (max)
 - Convex: 177°
 - Linear: 6.48 cm
 - Extended FOV
 - Convex: 217°
 - Linear: 40°
 - Convex Radius: 10 mm (Convex)
 - Steered Angle
 - B: -12°, -6°, 0°, 6°, 12°
 - C/PW: -15°-15°
 - Depth
 - Convex: 1.5-28.0 cm
 - Linear: 1.5-35.0 cm
 - Physical Footprint
 - Convex: 20 mm × 20 mm
 - Linear: 19.5 mm × 19.5 mm
 - Footprint
 - Convex: 20mm × 9mm

- Linear: 71.4mm × 9mm
- B-mode Frequencies
 - Convex: 3.5-6.0, 4.5-7.5, 5.5-8.8, 6.5-9.5 MHz
 - Linear: 4.8-7.8, 5.3-9.0, 5.8-10.0, 7.0-11.0 MHz
- Harmonic Frequencies
 - Convex: 8.0, 9.0, 10.0 MHz
 - Linear: 10.0, 11.0, 12.0 MHz
- Color Frequencies
 - Convex: 4.4, 5.0, 5.2, 5.5 (HR Flow) MHz
 - Linear: 6.0, 7.2, 8.3, 9.0 (HR Flow) MHz
- PW Frequencies
 - Convex: 4.5, 5.0, 5.5 MHz
 - Linear: 5.0, 6.3, 7.1, 8.3 MHz
- Biopsy Guide: NGB-051, multi angle, reusable; CIVCO 658-007 EX3, Biopsy Grid

10.6 Laparoscopic

LAP13-4Cs

- Application: Laparoscopic
- Advanced Function: Free Xros M, Color M, Contrast Imaging
- Bandwidth: 3.0-11.0 MHz
- Number of Elements: 128
- FOV (max): 30°
- Extended FOV: 70°
- Convex Radius: 50 mm
- Depth: 1.5-28.0 cm
- Physical Footprint: 10 mm × 10 mm
- Footprint: 31 mm × 6.5 mm
- B-mode Frequencies: 3.0-7.0, 4.0-9.0, 5.0-11.0 MHz
- Harmonic Frequencies: 8.0, 9.0, 10.0 MHz
- Color Frequencies: 4.4, 5.0, 5.7, 5.5 (HR Flow) MHz
- PW Frequencies: 4.4, 5.0, 5.7 MHz
- Biopsy Guide: not available

10.7 Pencil transducer

- CW2s
 - Application: Cardiac, Cephalic, Pediatric
 - Advanced Function: CW
 - Number of Elements: 2

- Biopsy Guide: not available
- CW5s
 - Application: Vascular, Cephalic, Pediatric
 - Advanced Function: CW
 - Number of Elements: 2
 - Biopsy Guide: not available

11 Peripheral Devices and Accessories

- Black/White Video Printer
 - Digital: MITSUBISHI P95DW-N
 - Digital& Analog: SONY UP-X898MD
- Color Digital Printer
SONY UP-D25MD
- Graph/Text Printer
HP OFFICEJET PRO 8100
Epson L130
Epson L805
Epson L3256
Epson L8058
Epson L3218
HP M283
HP 203dn
HP M208DW
HP P1108
HP 181 fw
HP M404DN
HP 7000(C9299A)
HP DeskJet 2029
HP LaserJet 1020 plus
HP Color Laserjet Pro M154a
HP DeskJet 1112
Canon G580
Canon G1830
Canon G2830
Canon TS708
Lexmark 2236
PANTUM P2550
- Gel Warmer
 - Support gel warming with 3 angle position: 15, 45, 90 degrees
 - Easily be disassembled off system for cleaning
 - Temperature with 4 levels: off/ 34°C/37°C/40°C, with deviation of ±2°C

- Light indicator for temperature protecting
- Dimension: 78mm (Width)*82mm (Depth)*119mm (Height)
- Weight: approx. 240g (net)
- Continuous operation time: >12H
- Footswitch
 - USB port: FS-81-SP-2(single pedal), 971-SWNOM (2/3 pedal)
 - Support user-definable functions (Freeze, Save, Print)
- ECG
 - 6-pin, AHA/IEC, for 3-lead wires
 - ECG wave display: on/off
 - ECG source: Lead/External
 - Position: 0 - 100%, 5%/step
 - Trig mode: off/single/dual/timer
 - Gain: 0-30, 1/step
 - Sweep speed: 1-6, 1/step
 - Invert: on/off
- PCG (not for sale in EU countries)
 - PCG wave display: on/off
 - Gain: 0-30, 1/step
 - Speed: 6 steps
 - Smooth: 1-4, 1/step
- Barcode Reader
 - SYMBOL LS2208 (1D)
 - SYMBOL DS4308 (2D)
- Built-in Wireless Adapter
 - Encryption: WPA, WPA2
 - Protocols: IEEE 802.11 ac/a/b/g/n
 - Frequency: 2.4G/5G
- iVocal Microphone
SAMSON XPD2

12 System Inputs and Outputs

- Audio input/output
Microphone: 1 port
- Video Output
 - S-Video out: 1 port
 - HDMI: 1 Port
 - VGA out: 1 port
- Physio Input
 - Support ECG/PCG signal
 - ECG: 1 port
 - PCG: 1 port
- Other Input/Output

- USB: 6 ports (5 USB 3.0 and 1 Type-C)
- Ethernet: 1 port

13 Safety and Conformance

- Quality Standards
 - ISO 9001
 - ISO 13485
- Design Standards
 - CSA C22.2 No. 601-1
 - EN 60601-1 and IEC 60601-1
 - EN 60601-1-2 and IEC 60601-1-2
 - EN 60601-1-6 and IEC 60601-1-6
 - EN 60601-2-37 and IEC 60601-2-37
 - EN 62304 and IEC 62304
 - EN 62366 and IEC 62366
 - EN ISO 17664 and ISO 17664

14 CE Declaration

The ultrasound system is fully in conformance with the Council Directive 93/42/EEC Concerning Medical Devices. The number adjacent to the CE marking (0123) is the code of the EU-notified body that certified meeting the requirements of Annex II excluding (4). of the Directive.

15 NOTICE

Not all features or specifications described in this document may be available in all transducers and/or modes.

Mindray reserves the right to make changes in specifications and features shown herein, or discontinue the product at any time without notice or obligation. Contact Mindray Representative for the most current information.

16 Appendix

16.1 Generic measurements

- Summarized items:
 - B-Mode: Distance, Angle, Circumference, Area, Volume, etc.
 - M-Mode: Slope, Heart Rate, Distance, Time, Velocity etc.

- D-Mode: PS/ED, Velocity, Heart Rate, Time, Acceleration, Velocity Ratio etc.
 - Detailed items:
 - B-Mode
 - Depth
 - Distance
 - Distance L-L
 - Shuttle
 - Ellipse
 - Trace
 - Spline
 - Cross
 - Angle (2-Lines)
 - Angle (3-Points)
 - Double Dist
 - Trace Len
 - Trace Len(Spline)
 - Parallel
 - Distance P-L
 - IMT
 - B-Profile
 - B-Hist(Ellipse)
 - B-Hist(Trace)
 - B-Hist(Spline)
 - B-Hist(Rectangle)
 - Color Vel
 - Strain Hist
 - Elas. Hist
 - Color Vel Profile
 - Elas.
 - Strain
 - TSM
 - Color Pixel Percent(Ellipse)
 - Color Pixel Percent(Trace)
 - Color Pixel Percent(Rectangle)
 - Color Pixel Percent(Recall)
 - Smart Calc
 - Smart Trace
 - Att.
 - Att. Hist
 -
 - Volume
 - Volume (Ellipse)
 - Volume (E+Dist.)
 - Ratio(D)
 - B Ratio
 -
 - M-Mode
 - HR
 - HR(R-R)
 - Slope
 - Distance
 - Depth
 - Time
 - Velocity
 - D-Mode
 - PS
 - ED
 - PS/ED
 - Vel
 - HR
 - HR(R-R)
 - Time
 - Auto Trace
 - Manual Trace
 - Spline Trace
 - Acceleration
- Volume
 - Volume
 - Volume (Ellipse)
 - Volume (E+Dist.)
 - Ratio(A)
 - Area1
 - Area2
 - Directional Ratio
 - D1
 - D2
 - RAC
 - Sag
 - XS
 - Volume Flow
 - Vas Area
 - TAMEAN
 - TAMAX
 - Elas. Ratio
 - A
 - B
 - Strain Ratio
 - A
 - B
 - Att. Ratio
 - Att A
 - Att B

 Ratio(Vel)
 Ratio(VTI)

 Volume Flow
 Vas Area
 TAMEAN
 TAMAX

16.2 Clinical option measurement package

- Abdomen
 Summarized items:
 Provide measurements of abdominal artery, abdominal vein, hepatic vein, liver, gallbladder, pancreas, appendix, pylorus, kidney, ureter, bladder, spleen, inferior vena cava, etc.

Detailed items:

B-Mode
 Aorta Bif
 AAA Status
 Shunt Diam
 Portal V Diam
 M Portal V Diam
 Splenic V Diam
 PS Conflnc Diam
 Renal V Diam
 SMV Diam
 IMV Diam
 CHD
 GB L
 GB H
 GB W
 GB wall th
 Cystic Duct
 CBD
 Panc duct
 Panc head
 Panc neck
 Panc body
 Panc tail
 Appendix
 Appendix Wall
 Pylorus
 Pylorus Wall
 Renal L
 Renal H

Renal W
 Cortex
 Adrenal L
 Adrenal H
 Adrenal W
 Ureter
 Cortex(Renal Transplant1)
 Renal V Diam(Renal Transplant1)
 Ureter Diam(Renal Transplant1)
 Cortex(Renal Transplant2)
 Renal V Diam(Renal Transplant2)
 Ureter Diam(Renal Transplant2)
 Smart Bladder Dual
 Pre-BL L
 Pre-BL H
 Pre-BL W
 Post-BL L
 Post-BL H
 Post-BL W
 Spleen L
 Spleen H
 Spleen W
 Spleen Area
 Skin-L.Capsule Dist.
 Hepatic Lesion1 Elas.
 Hepatic Lesion2 Elas.
 Hepatic Lesion3 Elas.
 LSM
 Rt DT(Insp)
 Rt DT(Expir)
 Lt DT(Insp)
 Lt DT(Expir)
 Free Fluid
 Smart HRI

 Renal Vol
 Pre-BL Vol
 Post-BL Vol
 Mictur.Vol

 Aorta
 Anterior-Posterior
 Transverse
 Outer Diameter
 Inner Diameter
 Outer Area
 Inner Area
 Celiac Axis

Anterior-Posterior	Long
Transverse	Anterior-Posterior
SMA	Transverse
Anterior-Posterior	GDA Aneurysm
Transverse	Long
C Hepatic A	Anterior-Posterior
Anterior-Posterior	Transverse
Transverse	IMA Aneurysm
Proper Hepatic A	Long
Anterior-Posterior	Anterior-Posterior
Transverse	Transverse
Hepatic A	EVAR Residual Aneurysm Sac(2D)
Anterior-Posterior	Anterior-Posterior
Transverse	Transverse
Splenic A	EVAR Inflow(2D)
Anterior-Posterior	Anterior-Posterior
Transverse	Transverse
GDA	EVAR Graft Body(2D)
Anterior-Posterior	Anterior-Posterior
Transverse	Transverse
IMA	EVAR Limb(2D)
Anterior-Posterior	Anterior-Posterior
Transverse	Transverse
AAA	EVAR Outflow(2D)
Long	Anterior-Posterior
Anterior-Posterior	Transverse
Transverse	Aortic Bypass Graft Anast(2D)
Celiac A Aneurysm	Anterior-Posterior
Long	Transverse
Anterior-Posterior	Aortic Bypass Graft Graft(2D)
Transverse	Anterior-Posterior
SMA Aneurysm	Transverse
Long	ABD Stenosis 1(2D)
Anterior-Posterior	Anterior-Posterior
Transverse	Transverse
C Hepatic A Aneurysm	Outer Diameter
Long	Inner Diameter
Anterior-Posterior	Outer Area
Transverse	Inner Area
Proper Hepatic A Aneurysm	ABD Stenosis 2(2D)
Long	Anterior-Posterior
Anterior-Posterior	Transverse
Transverse	Outer Diameter
Hepatic A Aneurysm	Inner Diameter
Long	Outer Area
Anterior-Posterior	Inner Area
Transverse	ABD Stenosis 3(2D)
Splenic A Aneurysm	Anterior-Posterior

Transverse	Hepatic Lesion 3
Outer Diameter	d1
Inner Diameter	d2
Outer Area	d3
Inner Area	Hepatic Cyst 1
ABD Stenosis 4(2D)	d1
Anterior-Posterior	d2
Transverse	d3
Outer Diameter	Hepatic Cyst 2
Inner Diameter	d1
Outer Area	d2
Inner Area	d3
IVC	Hepatic Cyst 3
Anterior-Posterior	d1
Transverse	d2
Checklist	d3
Hepatic V(2D)	GB
Anterior-Posterior	GB L
Transverse	GB H
Lt Hepatic V(2D)	GB W
Anterior-Posterior	GB wall th
Transverse	GB Finding 1
M Hepatic V(2D)	d1
Anterior-Posterior	d2
Transverse	d3
Rt Hepatic V(2D)	GB Finding 2
Anterior-Posterior	d1
Transverse	d2
Liver	d3
L	GB Finding 3
H	d1
W	d2
R Liver Lobe	d3
L	GB Finding 4
H	d1
W	d2
L Liver Lobe	d3
L	GB Finding 5
H	d1
W	d2
Hepatic Lesion 1	d3
d1	Panc Finding 1
d2	d1
d3	d2
Hepatic Lesion 2	d3
d1	Panc Finding 2
d2	d1
d3	d2

d3	W
Panc Finding 3	Kidney(Mid)
d1	H
d2	W
d3	Kidney(Inferior)
Panc Finding 4	H
d1	W
d2	Renal A
d3	Long
Panc Finding 5	Anterior-Posterior
d1	Transverse
d2	Renal A Aneurysm
d3	Long
Kidney	Anterior-Posterior
Renal L	Transverse
Renal H	Kidney(Renal Transplant1)
Renal W	L
Cortex	H
Adrenal	W
Adrenal L	Adrenal(Renal Transplant1)
Adrenal H	L
Adrenal W	H
Renal Lesion 1	W
d1	Finding 1(Renal Transplant1)
d2	L
d3	H
Renal Lesion 2	W
d1	Finding 2(Renal Transplant1)
d2	L
d3	H
Renal Lesion 3	W
d1	Finding 3(Renal Transplant1)
d2	L
d3	H
Renal Cyst 1	W
d1	Finding 4(Renal Transplant1)
d2	L
d3	H
Renal Cyst 2	W
d1	Finding 5(Renal Transplant1)
d2	L
d3	H
Renal Cyst 3	W
d1	Finding 6(Renal Transplant1)
d2	L
d3	H
Kidney(Superior)	W
H	Renal Transplant 1(2D)

Cortex(Renal Transplant1)	L
Renal V Diam(Renal Transplant1)	H
Ureter Diam(Renal Transplant1)	W
Kidney(Renal Transplant1)	Finding 3(Renal Transplant2)
L	L
H	H
W	W
Adrenal(Renal Transplant1)	Finding 4(Renal Transplant2)
L	L
H	H
W	W
Finding 1(Renal Transplant1)	Finding 5(Renal Transplant2)
L	L
H	H
W	W
Finding 2(Renal Transplant1)	Finding 6(Renal Transplant2)
L	L
H	H
W	W
Finding 3(Renal Transplant1)	Renal Transplant 2(2D)
L	Cortex(Renal Transplant2)
H	Renal V Diam(Renal Transplant2)
W	Ureter Diam(Renal Transplant2)
Finding 4(Renal Transplant1)	Kidney(Renal Transplant2)
L	L
H	H
W	W
Finding 5(Renal Transplant1)	Adrenal(Renal Transplant2)
L	L
H	H
W	W
Finding 6(Renal Transplant1)	Finding 1(Renal Transplant2)
L	L
H	H
W	W
Kidney(Renal Transplant2)	Finding 2(Renal Transplant2)
L	L
H	H
W	W
Adrenal(Renal Transplant2)	Finding 3(Renal Transplant2)
L	L
H	H
W	W
Finding 1(Renal Transplant2)	Finding 4(Renal Transplant2)
L	L
H	H
W	W
Finding 2(Renal Transplant2)	Finding 5(Renal Transplant2)

L	Hepatic A
H	Splenic A
W	GDA
Finding 6(Renal Transplant2)	IMA
L	Aorta(Post)
H	Celiac Axis(Post)
W	SMA(Post)
Bladder	C Hepatic A(Post)
Pre-BL L	Proper Hepatic A(Post)
Pre-BL H	Hepatic A(Post)
Pre-BL W	Splenic Artery(Post)
Post-BL L	GDA(Post)
Post-BL H	IMA(Post)
Post-BL W	EVAR Residual Aneurysm Sac
Smart Bladder	EVAR Inflow
XS Bladder	EVAR Graft Body
Sag Bladder	EVAR Limb
Spleen	EVAR Outflow
Spleen L	Aortic Bypass Graft Anast
Spleen H	Aortic Bypass Graft Graft
Spleen W	IVC Reflux
Spleen Area	IVC
Hepatic Lesion1 ElasRatio	Hepatic V
A	Lt Hepatic V
B	M Hepatic V
Hepatic Lesion2 ElasRatio	Rt Hepatic V
A	Portal V
B	M Portal V
Hepatic Lesion3 ElasRatio	Splenic V
A	Renal V
B	SMV
M-Mode	IMV
Rt DT(Insp M)	Hepatic A Anast(Liver Transplant)
Rt DT(Expir M)	Hepatic V Anast(Liver Transplant)
Lt DT(Insp M)	Portal V Anast(Liver Transplant)
Lt DT(Expir M)	IVC(Liver Transplant)
RDE(QB)	Hep V Confl(Liver Transplant)
RDE(DB)	Donor IVC(Liver Transplant)
LDE(QB)	Renal A
LDE(DB)	Ren A Org
D-Mode	M Renal A
Aorta	Renal A1
Celiac Axis	Renal A2
SMA	Hilum
C Hepatic A	Interlobar A
Proper Hepatic A	Arcuate A
	Segment A
	Artery Anast(Renal Transplant1)

Artery Anast 2(Renal Transplant1)
 Vein Anast(Renal Transplant1)
 Vein Anast 2(Renal Transplant1)
 Renal A(Renal Transplant1)
 Renal A1(Renal Transplant1)
 Renal A2(Renal Transplant1)
 Hilum(Renal Transplant1)
 Interlobar A(Renal Transplant1)
 Arcuate A(Renal Transplant1)
 Segmental A(Renal Transplant1)
 Renal Vein 1(Renal Transplant1)
 Renal Vein 2(Renal Transplant1)
 Artery Anast(Renal Transplant2)
 Artery Anast 2(Renal Transplant2)
 Vein Anast(Renal Transplant2)
 Vein Anast 2(Renal Transplant2)
 Renal A(Renal Transplant2)
 Renal A1(Renal Transplant2)
 Renal A2(Renal Transplant2)
 Hilum(Renal Transplant2)
 Interlobar A(Renal Transplant2)
 Arcuate A(Renal Transplant2)
 Segmental A(Renal Transplant2)
 Renal Vein 1(Renal Transplant2)
 Renal Vein 2(Renal Transplant2)
 TIPS

 SMA/Ao

CA/Ao

ABD Stenosis 1

Pre Sten

Sten

Post Sten

ABD Stenosis 2

Pre Sten

Sten

Post Sten

ABD Stenosis 3

Pre Sten

Sten

Post Sten

ABD Stenosis 4

Pre Sten

Sten

Post Sten

Renal Transplant 1(Doppler)

Artery Anast(Renal Transplant1)

Artery Anast 2(Renal Transplant1)
 Vein Anast(Renal Transplant1)
 Vein Anast 2(Renal Transplant1)
 Renal A(Renal Transplant1)
 Renal A1(Renal Transplant1)
 Renal A2(Renal Transplant1)
 Hilum(Renal Transplant1)
 Interlobar A(Renal Transplant1)
 Arcuate A(Renal Transplant1)
 Segmental A(Renal Transplant1)
 Renal Vein 1(Renal Transplant1)
 Renal Vein 2(Renal Transplant1)

Renal Transplant 2(Doppler)

Artery Anast(Renal Transplant2)
 Artery Anast 2(Renal Transplant2)
 Vein Anast(Renal Transplant2)
 Vein Anast 2(Renal Transplant2)
 Renal A(Renal Transplant2)
 Renal A1(Renal Transplant2)
 Renal A2(Renal Transplant2)
 Hilum(Renal Transplant2)
 Interlobar A(Renal Transplant2)
 Arcuate A(Renal Transplant2)
 Segmental A(Renal Transplant2)
 Renal Vein 1(Renal Transplant2)
 Renal Vein 2(Renal Transplant2)

- **Cardiology**

Summarized items:

- Left ventricular function measurement:
- Left ventricular function measurement method: Simpson, Mod.Simpson, S-P Ellipse, B-P Ellipse, Bullet, Teichholz, Cube, Gibson
- LV Mass (Cube-2D), (A-L), (T-E)
- LA Vol (Simp), (A-L)
- MVA VTI, AVA (VTI)
- Cardiac output: LVOT, RVOT, MV, TV
- Quantitative Analysis of PISA MR, PISA AR, PISA TR, and PISA PR
- Qp/Qs
- LV TEI Index
- RV TEI Index
- Z score (3 years and younger)
- Z score (Under 18)

Detailed items:

B-Mode
RVAWd(2D)
RVAWs(2D)
RVDd(2D)
RVDs(2D)
IVSd(2D)
IVSs(2D)
LVIDd(2D)
LVIDs(2D)
LVPWd(2D)
LVPWs(2D)
Diastole(2D)
Systole(2D)
LVLd apical
LVLs apical
LVAd apical
LVAs apical
LVAd sax MV
LVAs sax MV
LVAd sax Endo
LVAd sax Epi
LV Major
LV Minor
LV Area(d)
LV Area(s)
HR(2D)
RA Major
RA Minor
RA Area
RA Vol(A4C)
RAP
RV Area(d)
RV Area(s)
RV Major
RV Minor
LA Diam(2D)
LA Major
LA Minor
LA Area
LVOT Diam
Ao Diam(2D)
ACS(2D)
AV Diam
Ao Isthmus(2D)
Ao Sinus Diam(2D)
Ao st junct(2D)
AVA

Ao Arch Diam(2D)
Ao Asc Diam(2D)
Ao Desc Diam(2D)
Duct Art Diam
Post Ductal
Pre Ductal
MCS(2D)
MV Diam
MV EPSS(2D)
MVA
TV Diam
TVA
PV Diam
RVOT Diam
MPA Diam(2D)
RPA Diam(2D)
LPA Diam(2D)
IVC Diam(Expir)
IVC Diam(Insp)
SVC Diam(Expir)
SVC Diam(Insp)
LCA Diam
RCA Diam
PEd(2D)
PEs(2D)
VSD Diam
ASD Diam
PDA Diam
PFO Diam
AutoEF
Rt DT(Insp)
Rt DT(Expir)
Lt DT(Insp)
Lt DT(Expir)

LA/Ao(2D)

LV(2D)
Diastole(2D)
Systole(2D)
IVSd(2D)
LVIDd(2D)
LVPWd(2D)
IVSs(2D)
LVIDs(2D)
LVPWs(2D)
HR(2D)
Simpson

A4Cd	LA Vol(A4C)
A4Cs	LA Vol(A-L)
A2Cd	LA apical
A2Cs	LAA(A2C)
HR(2D)	LAA(A4C)
Mod. Simpson	MVA(VTI)
LVLd apical	LVOT Diam
LVLs apical	LVOT VTI
LVAd sax MV	MV VTI
LVAs sax MV	AVA(VTI)
LVAd sax PM	LVOT Diam
LVAs sax PM	LVOT VTI
HR(2D)	AV VTI
S-P Ellipse	CO(LVOT)
LVLd apical	LVOT Diam
LVAd apical	LVOT VTI
LVLs apical	AV HR
LVAs apical	CO(RVOT)
HR(2D)	RVOT Diam
B-P Ellipse	RVOT VTI
LVIDd(2D)	PV HR
LVAd sax MV	CO(MV)
LVIDs(2D)	MV Diam
LVAs sax MV	MV VTI
LVAd apical	MV HR
LVAs apical	CO(TV)
HR(2D)	TV Diam
Bullet	TV VTI
LVLd apical	TV HR
LVLs apical	PISA MR
LVAd sax MV	MR Rad
LVAs sax MV	MR Als. Vel
HR(2D)	MR VTI
LV Mass(Cube-2D)	PISA AR
IVSd(2D)	AR Rad
LVIDd(2D)	AR Als. Vel
LVPWd(2D)	AR VTI
LV Mass(A-L)	PISA TR
LVLd apical	TR Rad
LVAd sax Epi	TR Als. Vel
LVAd sax Endo	TR VTI
LV Mass(T-E)	PISA PR
LVAd sax Epi	PR Rad
LVAd sax Endo	PR Als. Vel
a	PR VTI
d	Qp/Qs
LA Vol(Simp)	LVOT Diam
LA Vol(A2C)	LVOT VTI

RVOT Diam	Ao Diam(M)
RVOT VTI	Ao Isthmus(M)
Z-Scores (= 3Y) (2D)	Ao Sinus Diam(M)
AV Diam	Ao st junct(M)
Ao Sinus Diam	ACS(M)
Ao st junct	HR(M)
PV Diam	IVSd(M)
Ao Arch IA-LCA	IVSs(M)
Ao Arch LCA-LSA	LA Diam(M)
Ao Arch after LSA	LPA Diam(M)
Ao Isthmus	Diastole(M)
Thoracic Ao Diam	Systole(M)
IVC Diam	LVET(M)
MV Diam	LVIDd(M)
TV Diam	LVIDs(M)
MPA Diam	LVOT Diam
RPA Diam	LVPEP(M)
LPA Diam	LVPWd(M)
Z-Scores (<18Y) (2D)	LVPWs(M)
LV Area(d) A4C	MCS(M)
LV Area(s) A4C	MPA Diam(M)
LVIDd A4C(2D)	MV A Amp
LVIDs A4C(2D)	MV E Amp
LA AP Diam A4C	MV D-E Slope
LA LL Diam A4C	MV D-E Amp
LA Area A4C	MV E-F Slope
RA AP Diam A4C	MV EPSS(M)
RA LL Diam A4C	PEd(M)
RA Area A4C	PEs(M)
RV Area(d) A4C	RPA Diam(M)
RV Area(s) A4C	RVET(M)
RVd Major A4C	RVOT Diam
RVs Major A4C	RVPEP(M)
RVd Minor (basal) A4C	MAPSE
RVd Minor (midcavity) A4C	TAPSE
LV Area(d) A2C	MV ALL
LV Area(s) A2C	IVC Diam(Insp)(M)
LVIDd A2C(2D)	IVC Diam(Expir)(M)
LVIDs A2C(2D)	SVC Diam(Insp)(M)
	SVC Diam(Expir)(M)
M-Mode	Rt DT(Insp M)
RVAWd(M)	Rt DT(Expir M)
RVAWs(M)	Lt DT(Insp M)
RVDd(M)	Lt DT(Expir M)
RVDs(M)	RDE(QB)
Ao Arch Diam(M)	RDE(DB)
Ao Asc Diam(M)	LDE(QB)
Ao Desc Diam(M)	LDE(DB)

LA/Ao(M)

LV(M)

Diastole(M)

Systole(M)

IVSd(M)

LVIDd(M)

LVPWd(M)

IVSs(M)

LVIDs(M)

LVPWs(M)

HR(M)

LV Mass(Cube-M)

IVSd(M)

LVIDd(M)

LVPWd(M)

LV Tei Index(M)

MV C-O dur(M)

LVET(M)

Z-Scores (= 3Y) (M)

IVSd(M)

LVPWd(M)

Z-Scores (<18Y) (M)

LVIDd(M)

LVIDs(M)

D-Mode

MV Aa(lateral)

MV Aa(medial)

AAo Vmax

AV VTI

AV HR

AV Vmax

AR DecT

AR Time

AR PHT

AR Ved

AR Vmax

AR VTI

MV ARa(lateral)

MV ARa(medial)

ASD Vmax

AV AccT

AV DecT

Coarc Post-Duct

Coarc Pre-Duct

DAo Vmax

MV DRa(lateral)

MV DRa(medial)

MV Ea(lateral)

MV Ea(medial)

IVC Vel(Expir)

IVC Vel(Insp)

IVCT

LPA Vmax

LVET(Doppler)

LVOT AccT

LVOT VTI

LVOT Vmax

LVPEP(Doppler)

MPA Vmax

dP/dt

Tau(BAI)

MR VTI

MR Vmax

MS Vmax

MV A Dur

MV A Vel

MV A VTI

MV AccT

MV DecT

MV E Dur

MV E Vel

MV E VTI

IVRT

MV VTI

MV HR

MV Vmax

PVein A Dur

PVein A Vel

PVein D Vel

PVein D VTI

PVein DecT

PVein S Vel

PVein S VTI

PDA Vel(d)

PDA Vel(s)

PR PHT

PR VTI

PR Ved

PR Vmax

PR DecT

PV AccT

PV VTI

PV HR

PV Vmax
 RAP
 RPA Vmax
 RVET(Doppler)
 RVOT Vmax
 RVOT VTI
 RVPEP(Doppler)
 MV Sa(lateral)
 MV Sa(medial)
 SVC Vel(Expir)
 SVC Vel(Insp)
 TR VTI
 TR Vmax
 TV A Dur
 TV A Vel
 TV AccT
 TV DecT
 TV E Vel
 TV VTI
 TV HR
 TV Vmax
 VSD Vmax
 Hepatic V S Vel
 Hepatic V D Vel

 MV E/A
 MVA(PHT)
 TV E/A
 TVA(PHT)

 LV Tei Index(Doppler)
 MV C-O dur(Doppler)
 LVET(Doppler)
 RVSP
 TR Vmax
 RAP
 PAEDP
 PR Ved
 RAP
 MVA(VTI)
 LVOT Diam
 LVOT VTI
 MV VTI
 AVA(VTI)
 LVOT Diam
 LVOT VTI
 AV VTI
 CO(LVOT)

LVOT Diam
 LVOT VTI
 AV HR
 CO(RVOT)
 RVOT Diam
 RVOT VTI
 PV HR
 CO(MV)
 MV Diam
 MV VTI
 MV HR
 CO(TV)
 TV Diam
 TV VTI
 TV HR
 RV Tei Index
 TV C-O dur
 RVET(Doppler)
 PISA MR
 MR Rad
 MR Als. Vel
 MR VTI
 PISA AR
 AR Rad
 AR Als. Vel
 AR VTI
 PISA TR
 TR Rad
 TR Als. Vel
 TR VTI
 PISA PR
 PR Rad
 PR Als. Vel
 PR VTI
 Qp/Qs
 LVOT Diam
 LVOT VTI
 RVOT Diam
 RVOT VTI
 • Emergency
 Summarized items: Provides
 application measurement related to
 emergency

 Detailed items:
 Renal L
 Renal H
 Renal W

CBD
 Portal V Diam
 CHD
 GB wall th
 Aorta Bif
 Ureter
 Pre-BL L
 Pre-BL H
 Pre-BL W
 Post-BL L
 Post-BL H
 Post-BL W
 Smart Bladder Dual
 GS
 YS L
 CRL
 BPD
 UT L
 UT H
 UT W
 Endo
 Ovary L
 Ovary H
 Ovary W

 Renal Vol
 Pre-BL Vol
 Post-BL Vol
 Mictur.Vol
 Ovary Vol
 UT Vol
 UT SUM

 Smart Bladder
 XS Bladder
 Sag Bladder
 Uterus
 UT L
 UT H
 UT W
 Endo
 Ovary
 Ovary L
 Ovary H
 Ovary W
 Kidney
 Renal L
 Renal H

Renal W
 Cortex
 Bladder
 Pre-BL L
 Pre-BL H
 Pre-BL W
 Post-BL L
 Post-BL H
 Post-BL W

M-Mode
 FHR (M)

D-Mode
 FHR (Doppler)

- Gynecology
 Summarized items:
 Provides measurements of uterus, cervix, ovary, follicle, urethra, rectum, levator ani muscle and sphincter.

Detailed items:

B-Mode
 UT L
 UT H
 UT W
 Endo
 Cervix L
 Cervix H
 Cervix W
 Ovary L
 Ovary H
 Ovary W
 Follicle1 L
 Follicle1 W
 Follicle1 H
 Follicle2 L
 Follicle2 W
 Follicle2 H
 Follicle3 L
 Follicle3 W
 Follicle3 H
 Follicle4 L
 Follicle4 W
 Follicle4 H
 Follicle5 L
 Follicle5 W
 Follicle5 H

Follicle6 L
Follicle6 W
Follicle6 H
Follicle7 L
Follicle7 W
Follicle7 H
Follicle8 L
Follicle8 W
Follicle8 H
Follicle9 L
Follicle9 W
Follicle9 H
Follicle10 L
Follicle10 W
Follicle10 H
Follicle11 L
Follicle11 W
Follicle11 H
Follicle12 L
Follicle12 W
Follicle12 H
Follicle13 L
Follicle13 W
Follicle13 H
Follicle14 L
Follicle14 W
Follicle14 H
Follicle15 L
Follicle15 W
Follicle15 H
Follicle16 L
Follicle16 W
Follicle16 H
Follicle17 L
Follicle17 W
Follicle17 H
Follicle18 L
Follicle18 W
Follicle18 H
Follicle19 L
Follicle19 W
Follicle19 H
Follicle20 L
Follicle20 W
Follicle20 H
Follicle21 L
Follicle21 W
Follicle21 H

Follicle22 L
Follicle22 W
Follicle22 H
Follicle23 L
Follicle23 W
Follicle23 H
Follicle24 L
Follicle24 W
Follicle24 H
Follicle25 L
Follicle25 W
Follicle25 H
Follicle26 L
Follicle26 W
Follicle26 H
Follicle27 L
Follicle27 W
Follicle27 H
Follicle28 L
Follicle28 W
Follicle28 H
Follicle29 L
Follicle29 W
Follicle29 H
Follicle30 L
Follicle30 W
Follicle30 H
DWT
BSD(R)
BSD(Va)
RVA(R)
RVA(Va)
UTA(R)
UTA(Va)
URA
PVA(R)
PVA(Va)
PUA(R)
PUA(Va)
BPW-SP Dist.(R)
BPW-SP Dist.(Va)
Cx-SP Dist.(R)
Cx-SP Dist.(Va)
RA-SP Dist.(R)
RA-SP Dist.(Va)
Shuttle(R)
Shuttle(Va)
Rectocele Depth

Intus. Depth	Follicle10
ARA(R)	Follicle11
ARA(Va)	Follicle12
ARA(C)	Follicle13
LH AP Diam(R)	Follicle14
LH AP Diam(Va)	Follicle15
LH AP Diam(C)	Follicle16
LH Lateral Diam(R)	Follicle17
LH Lateral Diam(Va)	Follicle18
LH Lateral Diam(C)	Follicle19
LH Area(R)	Follicle20
LH Area(Va)	Follicle21
LH Area(C)	Follicle22
LA Angle(R)	Follicle23
LA Angle(Va)	Follicle24
LA Angle(C)	Follicle25
LA Thickness(R)	Follicle26
LA Thickness(Va)	Follicle27
LA Thickness(C)	Follicle28
LUG(R)	Follicle29
LUG(Va)	Follicle30
LUG(C)	Mean DWT
GYN Lesion1 Strain	BND
GYN Lesion2 Strain	IAS Damage
GYN Lesion3 Strain	EAS Damage
Lesion1 Elas.	-----
Lesion2 Elas.	Uterus
Lesion3 Elas.	UT L
Fibroid1 Strain	UT H
Fibroid2 Strain	UT W
Fibroid3 Strain	Endo
Fibroid1 Elas.	Uterine Cervix
Fibroid2 Elas.	Cervix L
Fibroid3 Elas.	Cervix H
-----	Cervix W
UT Vol	Fibroid 1
UT SUM	d1
UT-L/CX-L	d2
Ovary Vol	d3
Follicle1	Fibroid 2
Follicle2	d1
Follicle3	d2
Follicle4	d3
Follicle5	Fibroid 3
Follicle6	d1
Follicle7	d2
Follicle8	d3
Follicle9	Uterine Finding 1

d1	d1
d2	d2
d3	d3
Uterine Finding 2	Ovarian Finding 4
d1	d1
d2	d2
d3	d3
Uterine Finding 3	Ovarian Finding 5
d1	d1
d2	d2
d3	d3
Uterine Finding 4	Ovarian Finding 6
d1	d1
d2	d2
d3	d3
Uterine Finding 5	Follicle1
d1	Follicle1 L
d2	Follicle1 W
d3	Follicle1 H
Uterine Finding 6	Follicle2
d1	Follicle2 L
d2	Follicle2 W
d3	Follicle2 H
Ovary	Follicle3
Ovary L	Follicle3 L
Ovary H	Follicle3 W
Ovary W	Follicle3 H
Ovarian Cyst 1	Follicle4
d1	Follicle4 L
d2	Follicle4 W
d3	Follicle4 H
Ovarian Cyst 2	Follicle5
d1	Follicle5 L
d2	Follicle5 W
d3	Follicle5 H
Ovarian Cyst 3	Follicle6
d1	Follicle6 L
d2	Follicle6 W
d3	Follicle6 H
Ovarian Finding 1	Follicle7
d1	Follicle7 L
d2	Follicle7 W
d3	Follicle7 H
Ovarian Finding 2	Follicle8
d1	Follicle8 L
d2	Follicle8 W
d3	Follicle8 H
Ovarian Finding 3	Follicle9

Follicle9 L	Follicle21 L
Follicle9 W	Follicle21 W
Follicle9 H	Follicle21 H
Follicle10	Follicle22
Follicle10 L	Follicle22 L
Follicle10 W	Follicle22 W
Follicle10 H	Follicle22 H
Follicle11	Follicle23
Follicle11 L	Follicle23 L
Follicle11 W	Follicle23 W
Follicle11 H	Follicle23 H
Follicle12	Follicle24
Follicle12 L	Follicle24 L
Follicle12 W	Follicle24 W
Follicle12 H	Follicle24 H
Follicle13	Follicle25
Follicle13 L	Follicle25 L
Follicle13 W	Follicle25 W
Follicle13 H	Follicle25 H
Follicle14	Follicle26
Follicle14 L	Follicle26 L
Follicle14 W	Follicle26 W
Follicle14 H	Follicle26 H
Follicle15	Follicle27
Follicle15 L	Follicle27 L
Follicle15 W	Follicle27 W
Follicle15 H	Follicle27 H
Follicle16	Follicle28
Follicle16 L	Follicle28 L
Follicle16 W	Follicle28 W
Follicle16 H	Follicle28 H
Follicle17	Follicle29
Follicle17 L	Follicle29 L
Follicle17 W	Follicle29 W
Follicle17 H	Follicle29 H
Follicle18	Follicle30
Follicle18 L	Follicle30 L
Follicle18 W	Follicle30 W
Follicle18 H	Follicle30 H
Follicle19	GYN Lesion 1
Follicle19 L	d1
Follicle19 W	d2
Follicle19 H	d3
Follicle20	GYN Lesion 2
Follicle20 L	d1
Follicle20 W	d2
Follicle20 H	d3
Follicle21	GYN Lesion 3

- d1
- d2
- d3
- Residual Urine
 - BL Height
 - BL Depth
- GYN Lesion1 Strain Ratio
 - A
 - B
- GYN Lesion2 Strain Ratio
 - A
 - B
- GYN Lesion3 Strain Ratio
 - A
 - B
- Lesion1 Elas. Ratio
 - A
 - B
- Lesion2 Elas. Ratio
 - A
 - B
- Lesion3 Elas. Ratio
 - A
 - B
- Fibroid1 Strain Ratio
 - A
 - B
- Fibroid2 Strain Ratio
 - A
 - B
- Fibroid3 Strain Ratio
 - A
 - B
- Fibroid1 Elas. Ratio
 - A
 - B
- Fibroid2 Elas. Ratio
 - A
 - B
- Fibroid3 Elas. Ratio
 - A
 - B
- Obstetric
 - Summarized items:
 - Multiple fetuses, maximum number of fetuses: 4
 - Fetal physiology score
 - Z score
 - Fetal GA
 - Fetal growth curve
 - EFW
 - Detailed items:
 - B-Mode
 - GS
 - YS L
 - CRL
 - NT
 - BPD
 - OFD
 - HC
 - AC
 - FL
 - TAD
 - APAD
 - TCD
 - CM
 - IT
 - LVW
 - HW
 - OOD
 - IOD
 - HUM
 - Ulna
 - RAD
 - Tibia
 - FIB
 - CLAV
 - Vertebrae
 - MP
 - Foot
 - NBL
 - Ear
 - APTD
 - TTD
 - FTA
 - THD
 - HrtC
 - TC
 - Umb VD
 - F-kidney L
 - Mat Kidney
 - Cervix L
 - AF
 - NF
 - Orbit
 - PL Thickness

Sac Diam1	Endo
Sac Diam2	AH
Sac Diam3	PH
AF1	3th Ventricle
AF2	NT Above Cord
AF3	NT Below Cord
AF4	Mandible
LVIDd	Prenasal th
LVIDs	Heart AP
LV Diam	Heart T
LA Diam	LV Width
RVIDd	LV Length
RVIDs	RV Width
RV Diam	RV Length
RA Diam	LA Width
IVSd	RA Width
IVSs	LVWd
IVS	LVWs
LV Area	RVWd
LA Area	RVWs
RV Area	AV Diam
RA Area	AV Area
Ao Diam	PV Area
MPA Diam	F-kidney H
LVOT Diam	F-kidney W
RVOT Diam	Lung
Facial Angle	Stomach
HrtA	YS H
MV Diam(Z-Score)	YS W
PV Diam(Z-Score)	Amniotic Sac L
Ao Asc Diam(Z-Score)	Amniotic Sac H
Ao Desc Diam(Z-Score)	Amniotic Sac W
Duct Art Diam(Z-Score)	Ovary Cyst L
TV Diam(Z-Score)	Ovary Cyst H
LPA Diam(Z-Score)	Ovary Cyst W
RPA Diam(Z-Score)	UT AW
IVC Diam(Z-Score)	UT PW
AV Diam(Z-Score)	CSP
MPA Diam(Z-Score)	FMF
RV Diam(Z-Score)	MMF
LV Diam(Z-Score)	Lung CCAM L
RV Area(Z-Score)	Lung CCAM H
LV Area(Z-Score)	Lung CCAM W
RVIDd(Z-Score)	AD
LVIDd(Z-Score)	Iliac Wing Angle
UT L	FAGL
UT H	FAG
UT W	Intestinum Crassum

Liver Length
 Rib Length
 Shoulder Blade

 MAD
 Mean Sac Diam
 AFI
 EFW
 EFW2
 HC/AC(Campbell)
 FL/AC
 FL/BPD
 AXT
 CI
 FL/HC(Hadlock)
 AC(c)
 HC(c)
 HrtC/TC
 TCD/AC
 LVW/HW
 LVD/RVD
 LAD/RAD
 AoD/MPAD
 LAD/AoD
 UT Vol
 UT SUM
 UT-L/CX-L

 AFI
 AF1
 AF2
 AF3
 AF4
 Uterus
 UT L
 UT H
 UT W
 Endo

 M-Mode
 FHR (M)
 LVIDd
 LVIDs
 RVIDd
 RVIDs
 IVSd
 IVSs
 RVIDd(Z-Score)

LVIDd(Z-Score)
 MVE
 TVE
 AVE
 MAPSE
 TAPSE
 LV ICT
 LV IRT
 LV ET
 RV ICT
 RV IRT
 RV ET

 D-Mode
 Umb A
 Duct Veno
 Placenta A
 MCA
 Fetal Ao
 Desc Aorta
 Ut A
 Ovarian A
 FHR (Doppler)
 Asc Aorta
 RVOT
 LVOT
 MV E
 MV A
 TV E
 TV A
 MV E'
 MV A'
 MV S'
 TV E'
 TV A'
 TV S'
 AV PV
 AV VTI
 PV PV
 PV VTI
 Duct Art PV
 Duct Art VTI
 AV TPV
 PV TPV
 Duct Art TPV
 Thoracic Aorta
 Hepatic Vein
 IVC

- Umb V
- Ovary
- Endometrium
- Cervical Cancer
- Fibroid
- Duct Art
- ICA
- Celiac A
-
- MV E/A
- TV E/A
- MV E/E'
- TV E/E'
- Pediatrics
 - Summarized items: provides application measurement related to pediatrics
 - Detailed items:
 - B-Mode
 - HIP
 - HIP-Graft
 - HIP(α)
 - HIP(β)
 - d/D
- SmallPart
 - Summarized items: Provides measurements of thyroid, parotid gland, lymph nodes, testis, epididymis, thyroid nodules, breast, etc.
 - Detailed items:
 - B-Mode
 - Thyroid L
 - Thyroid H
 - Thyroid W
 - Isthmus H
 - THY Nodule1 Strain
 - THY Nodule2 Strain
 - THY Nodule3 Strain
 - THY Nodule1 Elas.
 - THY Nodule2 Elas.
 - THY Nodule3 Elas.
 - Breast Mass1 Strain
 - Breast Mass1 Elas.
 - Breast Mass2 Strain
 - Breast Mass2 Elas.

- Breast Mass3 Strain
- Breast Mass3 Elas.
- Breast Mass4 Strain
- Breast Mass4 Elas.
- Breast Mass5 Strain
- Breast Mass5 Elas.
- Breast Mass6 Strain
- Breast Mass6 Elas.
- Breast Mass7 Strain
- Breast Mass7 Elas.
- Breast Mass8 Strain
- Breast Mass8 Elas.
- Breast Mass9 Strain
- Breast Mass9 Elas.
- Breast Mass10 Strain
- Breast Mass10 Elas.
- Testicular L
- Testicular H
- Testicular W
- Epididymis L
- Epididymis H
- Epididymis W
- Scrotal Wall
- Testis V(2D)
- Testis V(Valsalva 2D)
-
- Thyroid Vol
- Testicular Vol
-
- Thyroid(Superior)
 - Anterior-Posterior
 - Transverse
- Thyroid(Mid)
 - Anterior-Posterior
 - Transverse
- Thyroid(Inferior)
 - Anterior-Posterior
 - Transverse
- Parathyroid 1
 - Long
 - Anterior-Posterior
 - Transverse
- Parathyroid 2
 - Long
 - Anterior-Posterior
 - Transverse
- Parotid
 - Long

Anterior-Posterior
Transverse
Lymph Node 1
Long
Anterior-Posterior
Transverse
Lymph Node 2
Long
Anterior-Posterior
Transverse
Lymph Node 3
Long
Anterior-Posterior
Transverse
Lymph Node 4
Long
Anterior-Posterior
Transverse
Lymph Node 5
Long
Anterior-Posterior
Transverse
Lymph Node 6
Long
Anterior-Posterior
Transverse
Thyroid
Thyroid L
Thyroid H
Thyroid W
Thyroid Nodule 1
Long
Anterior-Posterior
Transverse
Thyroid Nodule 2
Long
Anterior-Posterior
Transverse
Thyroid Nodule 3
Long
Anterior-Posterior
Transverse
Thyroid Nodule 4
Long
Anterior-Posterior
Transverse
Thyroid Nodule 5
Long

Anterior-Posterior
Transverse
Thyroid Nodule 6
Long
Anterior-Posterior
Transverse
Thyroid Nodule 7
Long
Anterior-Posterior
Transverse
Thyroid Nodule 8
Long
Anterior-Posterior
Transverse
Thyroid Nodule 9
Long
Anterior-Posterior
Transverse
Thyroid Nodule 10
Long
Anterior-Posterior
Transverse
THY Nodule1 Strain Ratio
A
B
THY Nodule2 Strain Ratio
A
B
THY Nodule3 Strain Ratio
A
B
THY Nodule1 Elas. Ratio
A
B
THY Nodule2 Elas. Ratio
A
B
THY Nodule3 Elas. Ratio
A
B
Breast Mass 1
L
H
W
Nip. Dist.
Skin Dist.
Breast Mass 2
L

H
W
Nip. Dist.
Skin Dist.

Breast Mass 3

L

H

W

Nip. Dist.

Skin Dist.

Breast Mass 4

L

H

W

Nip. Dist.

Skin Dist.

Breast Mass 5

L

H

W

Nip. Dist.

Skin Dist.

Breast Mass 6

L

H

W

Nip. Dist.

Skin Dist.

Breast Mass 7

L

H

W

Nip. Dist.

Skin Dist.

Breast Mass 8

L

H

W

Nip. Dist.

Skin Dist.

Breast Mass 9

L

H

W

Nip. Dist.

Skin Dist.

Breast Mass 10

L

H

W

Nip. Dist.

Skin Dist.

Breast Mass1 Strain Ratio

A

B

Breast Mass1 Elas. Ratio

A

B

Breast Mass2 Strain Ratio

A

B

Breast Mass2 Elas. Ratio

A

B

Breast Mass3 Strain Ratio

A

B

Breast Mass3 Elas. Ratio

A

B

Breast Mass4 Strain Ratio

A

B

Breast Mass4 Elas. Ratio

A

B

Breast Mass5 Strain Ratio

A

B

Breast Mass5 Elas. Ratio

A

B

Breast Mass6 Strain Ratio

A

B

Breast Mass6 Elas. Ratio

A

B

Breast Mass7 Strain Ratio

A

B

Breast Mass7 Elas. Ratio

A

B

Breast Mass8 Strain Ratio

A

B
 Breast Mass8 Elas. Ratio
 A
 B
 Breast Mass9 Strain Ratio
 A
 B
 Breast Mass9 Elas. Ratio
 A
 B
 Breast Mass10 Strain Ratio
 A
 B
 Breast Mass10 Elas. Ratio
 A
 B
 Testicular
 Testicular L
 Testicular H
 Testicular W
 Testis Mass 1
 d1
 d2
 d3
 Testis Mass 2
 d1
 d2
 d3
 Testis Mass 3
 d1
 d2
 d3
 Epididymis
 Epididymis L
 Epididymis H
 Epididymis W
 Testicle(Superior)
 H
 W
 Testicle(Mid)
 H
 W
 Testicle(Inferior)
 H
 W
 Epididymal Head
 L
 H

W
 Epididymal Body
 L
 H
 W
 Epididymal Tail
 L
 H
 W

 M-Mode

 D-Mode
 STA
 ITA
 Isthmus
 Parathyroid 1
 Parathyroid 2
 Testis A
 Testis V
 Testis V(Valsalva)
 Epididymis A
 Epididymis V
 • Urology
 Summarized items:
 Provides measurement of kidney, kidney cortex, adrenal gland, prostate gland, seminal vesicle, bladder, testis, epididymis, scrotal wall, and urine output.

 Detailed items:
 B-Mode
 Renal L
 Renal H
 Renal W
 Cortex
 Adrenal L
 Adrenal H
 Adrenal W
 Ureter
 Cortex(Renal Transplant1)
 Renal V Diam(Renal Transplant1)
 Ureter Diam(Renal Transplant1)
 Cortex(Renal Transplant2)
 Renal V Diam(Renal Transplant2)
 Ureter Diam(Renal Transplant2)
 Prostate L

Prostate H	d3
Prostate W	Renal Lesion 2
Seminal L	d1
Seminal H	d2
Seminal W	d3
Urethra	Renal Lesion 3
Smart Bladder Dual	d1
Pre-BL L	d2
Pre-BL H	d3
Pre-BL W	Renal Cyst 1
Post-BL L	d1
Post-BL H	d2
Post-BL W	d3
Testicular L	Renal Cyst 2
Testicular H	d1
Testicular W	d2
Epididymis L	d3
Epididymis H	Renal Cyst 3
Epididymis W	d1
Scrotal Wall	d2
Testis V(2D)	d3
Testis V(Valsalva 2D)	Kidney(Superior)
Prostate Mass1 Strain	H
Prostate Mass2 Strain	W
Prostate Mass3 Strain	Kidney(Mid)
Prostate Mass1 Elas.	H
Prostate Mass2 Elas.	W
Prostate Mass3 Elas.	Kidney(Inferior)
-----	H
Renal Vol	W
Prostate Vol	Renal A
Pre-BL Vol	Long
Post-BL Vol	Anterior-Posterior
Mictur.Vol	Transverse
Testicular Vol	Renal A Aneurysm
-----	Long
Kidney	Anterior-Posterior
Renal L	Transverse
Renal H	Kidney(Renal Transplant1)
Renal W	L
Cortex	H
Adrenal	W
Adrenal L	Adrenal(Renal Transplant1)
Adrenal H	L
Adrenal W	H
Renal Lesion 1	W
d1	Finding 1(Renal Transplant1)
d2	L

H	H
W	W
Finding 2(Renal Transplant1)	Finding 5(Renal Transplant1)
L	L
H	H
W	W
Finding 3(Renal Transplant1)	Finding 6(Renal Transplant1)
L	L
H	H
W	W
Finding 4(Renal Transplant1)	Kidney(Renal Transplant2)
L	L
H	H
W	W
Finding 5(Renal Transplant1)	Adrenal(Renal Transplant2)
L	L
H	H
W	W
Finding 6(Renal Transplant1)	Finding 1(Renal Transplant2)
L	L
H	H
W	W
Renal Transplant 1(2D)	Finding 2(Renal Transplant2)
Cortex(Renal Transplant1)	L
Renal V Diam(Renal Transplant1)	H
Ureter Diam(Renal Transplant1)	W
Kidney(Renal Transplant1)	Finding 3(Renal Transplant2)
L	L
H	H
W	W
Adrenal(Renal Transplant1)	Finding 4(Renal Transplant2)
L	L
H	H
W	W
Finding 1(Renal Transplant1)	Finding 5(Renal Transplant2)
L	L
H	H
W	W
Finding 2(Renal Transplant1)	Finding 6(Renal Transplant2)
L	L
H	H
W	W
Finding 3(Renal Transplant1)	Renal Transplant 2(2D)
L	Cortex(Renal Transplant2)
H	Renal V Diam(Renal Transplant2)
W	Ureter Diam(Renal Transplant2)
Finding 4(Renal Transplant1)	Kidney(Renal Transplant2)
L	L

H	d2
W	d3
Adrenal(Renal Transplant2)	Prostate Mass 3
L	d1
H	d2
W	d3
Finding 1(Renal Transplant2)	Bladder
L	Pre-BL L
H	Pre-BL H
W	Pre-BL W
Finding 2(Renal Transplant2)	Post-BL L
L	Post-BL H
H	Post-BL W
W	Smart Bladder
Finding 3(Renal Transplant2)	XS Bladder
L	Sag Bladder
H	Testicular
W	Testicular L
Finding 4(Renal Transplant2)	Testicular H
L	Testicular W
H	Testis Mass 1
W	d1
Finding 5(Renal Transplant2)	d2
L	d3
H	Testis Mass 2
W	d1
Finding 6(Renal Transplant2)	d2
L	d3
H	Testis Mass 3
W	d1
Prostate	d2
Prostate L	d3
Prostate H	Epididymis
Prostate W	Epididymis L
Prostate2	Epididymis H
Long	Epididymis W
Anterior-Posterior	Testicle(Superior)
Coronal	H
Seminal Vesicle	W
Seminal L	Testicle(Mid)
Seminal H	H
Seminal W	W
Prostate Mass 1	Testicle(Inferior)
d1	H
d2	W
d3	Epididymal Head
Prostate Mass 2	L
d1	H

W	Interlobar A(Renal Transplant1)
Epididymal Body	Arcuate A(Renal Transplant1)
L	Segmental A(Renal Transplant1)
H	Renal Vein 1(Renal Transplant1)
W	Renal Vein 2(Renal Transplant1)
Epididymal Tail	Artery Anast(Renal Transplant2)
L	Artery Anast 2(Renal Transplant2)
H	Vein Anast(Renal Transplant2)
W	Vein Anast 2(Renal Transplant2)
Prostate Mass1 Strain Ratio	Renal A(Renal Transplant2)
A	Renal A1(Renal Transplant2)
B	Renal A2(Renal Transplant2)
Prostate Mass2 Strain Ratio	Hilum(Renal Transplant2)
A	Interlobar A(Renal Transplant2)
B	Arcuate A(Renal Transplant2)
Prostate Mass3 Strain Ratio	Segmental A(Renal Transplant2)
A	Renal Vein 1(Renal Transplant2)
B	Renal Vein 2(Renal Transplant2)
Prostate Mass1 Elas. Ratio	Testis A
A	Testis V
B	Testis V(Valsalva)
Prostate Mass2 Elas. Ratio	Epididymis A
A	Epididymis V
B	-----
Prostate Mass3 Elas. Ratio	Renal Transplant 1(Doppler)
A	Artery Anast(Renal Transplant1)
B	Artery Anast 2(Renal Transplant1)
	Vein Anast(Renal Transplant1)
	Vein Anast 2(Renal Transplant1)
	Renal A(Renal Transplant1)
	Renal A1(Renal Transplant1)
	Renal A2(Renal Transplant1)
	Hilum(Renal Transplant1)
	Interlobar A(Renal Transplant1)
	Arcuate A(Renal Transplant1)
	Segmental A(Renal Transplant1)
	Renal Vein 1(Renal Transplant1)
	Renal Vein 2(Renal Transplant1)
	Renal Transplant 2(Doppler)
	Artery Anast(Renal Transplant2)
	Artery Anast 2(Renal Transplant2)
	Vein Anast(Renal Transplant2)
	Vein Anast 2(Renal Transplant2)
	Renal A(Renal Transplant2)
	Renal A1(Renal Transplant2)
	Renal A2(Renal Transplant2)
	Hilum(Renal Transplant2)
	Interlobar A(Renal Transplant2)
M-Mode	
D-Mode	
Renal A	
Ren A Org	
M Renal A	
Renal A1	
Renal A2	
Hilum	
Interlobar A	
Arcuate A	
Segment A	
Artery Anast(Renal Transplant1)	
Artery Anast 2(Renal Transplant1)	
Vein Anast(Renal Transplant1)	
Vein Anast 2(Renal Transplant1)	
Renal A(Renal Transplant1)	
Renal A1(Renal Transplant1)	
Renal A2(Renal Transplant1)	
Hilum(Renal Transplant1)	

- Arcuate A(Renal Transplant2)
- Segmental A(Renal Transplant2)
- Renal Vein 1(Renal Transplant2)
- Renal Vein 2(Renal Transplant2)
- Vascular
 - Summarized items:
 - Provides measurement of carotid artery, jugular vein, upper limb artery, lower limb artery, upper limb vein and lower limb vein.
 - Detailed items:
 - B-Mode
 - CCA IMT
 - Bulb IMT
 - ICA IMT
 - ECA IMT
 -
 - IMT
 - CCA IMT
 - Bulb IMT
 - ICA IMT
 - ECA IMT
 - CCA
 - Anterior-Posterior
 - Transverse
 - Outer Diameter
 - Inner Diameter
 - Outer Area
 - Inner Area
 - Bulb
 - Anterior-Posterior
 - Transverse
 - Outer Diameter
 - Inner Diameter
 - Outer Area
 - Inner Area
 - Carotid Bifurcation
 - Anterior-Posterior
 - Transverse
 - Outer Diameter
 - Inner Diameter
 - Outer Area
 - Inner Area
 - ICA
 - Anterior-Posterior
 - Transverse
 - Outer Diameter

- Inner Diameter
- Outer Area
- Inner Area
- ECA
 - Anterior-Posterior
 - Transverse
 - Outer Diameter
 - Inner Diameter
 - Outer Area
 - Inner Area
- Vert A
 - Anterior-Posterior
 - Transverse
 - Outer Diameter
 - Inner Diameter
 - Outer Area
 - Inner Area
- Subclav A
 - Anterior-Posterior
 - Transverse
 - Outer Diameter
 - Inner Diameter
 - Outer Area
 - Inner Area
- Innom A
 - Anterior-Posterior
 - Transverse
 - Outer Diameter
 - Inner Diameter
 - Outer Area
 - Inner Area
- Mammary A
 - Anterior-Posterior
 - Transverse
 - Outer Diameter
 - Inner Diameter
 - Outer Area
 - Inner Area
- CCA Aneurysm
 - Long
 - Anterior-Posterior
 - Transverse
- Bulb Aneurysm
 - Long
 - Anterior-Posterior
 - Transverse
- Carotid Bifurcation Aneurysm
 - Long

Anterior-Posterior
 Transverse
 ICA Aneurysm
 Long
 Anterior-Posterior
 Transverse
 ECA Aneurysm
 Long
 Anterior-Posterior
 Transverse
 Vert A Aneurysm
 Long
 Anterior-Posterior
 Transverse
 Subclav A Aneurysm
 Long
 Anterior-Posterior
 Transverse
 Innom A Aneurysm
 Long
 Anterior-Posterior
 Transverse
 Mammary A Aneurysm
 Long
 Anterior-Posterior
 Transverse
 Carotid Graft 1 Anast
 Long
 Anterior-Posterior
 Transverse
 Carotid Graft 1 Graft
 Long
 Anterior-Posterior
 Transverse
 Carotid Graft 2 Anast
 Long
 Anterior-Posterior
 Transverse
 Carotid Graft 2 Graft
 Long
 Anterior-Posterior
 Transverse
 Carotid Graft 3 Anast
 Long
 Anterior-Posterior
 Transverse
 Carotid Graft 3 Graft
 Long

Anterior-Posterior
 Transverse
 Carotid Stent 1
 Long
 Anterior-Posterior
 Transverse
 Carotid Stent 2
 Long
 Anterior-Posterior
 Transverse
 Carotid Stent 3
 Long
 Anterior-Posterior
 Transverse
 Carotid Stenosis 1
 Anterior-Posterior
 Transverse
 Outer Diameter
 Inner Diameter
 Outer Area
 Inner Area
 Carotid Stenosis 2
 Anterior-Posterior
 Transverse
 Outer Diameter
 Inner Diameter
 Outer Area
 Inner Area
 Carotid Stenosis 3
 Anterior-Posterior
 Transverse
 Outer Diameter
 Inner Diameter
 Outer Area
 Inner Area
 Carotid Stenosis 4
 Anterior-Posterior
 Transverse
 Outer Diameter
 Inner Diameter
 Outer Area
 Inner Area
 Axill A
 Anterior-Posterior
 Transverse
 Outer Diameter
 Inner Diameter
 Outer Area

Inner Area	Anterior-Posterior
Brachial A	Transverse
Anterior-Posterior	UE A Graft 2 Native Inflow
Transverse	Anterior-Posterior
Outer Diameter	Transverse
Inner Diameter	UE A Graft 2 Anast
Outer Area	Anterior-Posterior
Inner Area	Transverse
Radial A	UE A Graft 2 Graft
Anterior-Posterior	Anterior-Posterior
Transverse	Transverse
Outer Diameter	UE A Graft 2 Native Outflow
Inner Diameter	Anterior-Posterior
Outer Area	Transverse
Inner Area	UE A Graft 3 Native Inflow
Ulnar A	Anterior-Posterior
Anterior-Posterior	Transverse
Transverse	UE A Graft 3 Anast
Outer Diameter	Anterior-Posterior
Inner Diameter	Transverse
Outer Area	UE A Graft 3 Graft
Inner Area	Anterior-Posterior
Axill A Aneurysm	Transverse
Long	UE A Graft 3 Native Outflow
Anterior-Posterior	Anterior-Posterior
Transverse	Transverse
Brachial A Aneurysm	UE A Stent 1
Long	Long
Anterior-Posterior	Anterior-Posterior
Transverse	Transverse
Radial A Aneurysm	UE A Stent 2
Long	Long
Anterior-Posterior	Anterior-Posterior
Transverse	Transverse
Ulnar A Aneurysm	UE A Stent 3
Long	Long
Anterior-Posterior	Anterior-Posterior
Transverse	Transverse
UE A Graft 1 Native Inflow	UE A Stenosis 1(2D)
Anterior-Posterior	Anterior-Posterior
Transverse	Transverse
UE A Graft 1 Anast	Outer Diameter
Anterior-Posterior	Inner Diameter
Transverse	Outer Area
UE A Graft 1 Graft	Inner Area
Anterior-Posterior	UE A Stenosis 2(2D)
Transverse	Anterior-Posterior
UE A Graft 1 Native Outflow	Transverse

Outer Diameter
 Inner Diameter
 Outer Area
 Inner Area
 UE A Stenosis 3(2D)
 Anterior-Posterior
 Transverse
 Outer Diameter
 Inner Diameter
 Outer Area
 Inner Area
 UE A Stenosis 4(2D)
 Anterior-Posterior
 Transverse
 Outer Diameter
 Inner Diameter
 Outer Area
 Inner Area
 C.Iliac A
 Anterior-Posterior
 Transverse
 Outer Diameter
 Inner Diameter
 Outer Area
 Inner Area
 Ex.Iliac A
 Anterior-Posterior
 Transverse
 Outer Diameter
 Inner Diameter
 Outer Area
 Inner Area
 IIA
 Anterior-Posterior
 Transverse
 Outer Diameter
 Inner Diameter
 Outer Area
 Inner Area
 CFA
 Anterior-Posterior
 Transverse
 Outer Diameter
 Inner Diameter
 Outer Area
 Inner Area
 DFA
 Anterior-Posterior

Transverse
 Outer Diameter
 Inner Diameter
 Outer Area
 Inner Area
 SFA
 Anterior-Posterior
 Transverse
 Outer Diameter
 Inner Diameter
 Outer Area
 Inner Area
 Pop A
 Anterior-Posterior
 Transverse
 Outer Diameter
 Inner Diameter
 Outer Area
 Inner Area
 TP Trunk A
 Anterior-Posterior
 Transverse
 Outer Diameter
 Inner Diameter
 Outer Area
 Inner Area
 A.Tib A
 Anterior-Posterior
 Transverse
 Outer Diameter
 Inner Diameter
 Outer Area
 Inner Area
 Peroneal A
 Anterior-Posterior
 Transverse
 Outer Diameter
 Inner Diameter
 Outer Area
 Inner Area
 P.Tib A
 Anterior-Posterior
 Transverse
 Outer Diameter
 Inner Diameter
 Outer Area
 Inner Area
 Dors.Ped. A

Anterior-Posterior
Transverse
Outer Diameter
Inner Diameter
Outer Area
Inner Area
C.Iliac A Aneurysm
Long
Anterior-Posterior
Transverse
Ex.Iliac A Aneurysm
Long
Anterior-Posterior
Transverse
IIA Aneurysm
Long
Anterior-Posterior
Transverse
CFA Aneurysm
Long
Anterior-Posterior
Transverse
DFA Aneurysm
Long
Anterior-Posterior
Transverse
SFA Aneurysm
Long
Anterior-Posterior
Transverse
Pop A Aneurysm
Long
Anterior-Posterior
Transverse
TP Trunk A Aneurysm
Long
Anterior-Posterior
Transverse
A.Tib A Aneurysm
Long
Anterior-Posterior
Transverse
Peroneal A Aneurysm
Long
Anterior-Posterior
Transverse
P.Tib A Aneurysm
Long

Anterior-Posterior
Transverse
Dors.Ped. A Aneurysm
Long
Anterior-Posterior
Transverse
LE A Graft 1 Native Inflow
Anterior-Posterior
Transverse
LE A Graft 1 Anast
Anterior-Posterior
Transverse
LE A Graft 1 Graft
Anterior-Posterior
Transverse
LE A Graft 1 Native Outflow
Anterior-Posterior
Transverse
LE A Graft 2 Native Inflow
Anterior-Posterior
Transverse
LE A Graft 2 Anast
Anterior-Posterior
Transverse
LE A Graft 2 Graft
Anterior-Posterior
Transverse
LE A Graft 2 Native Outflow
Anterior-Posterior
Transverse
LE A Graft 3 Native Inflow
Anterior-Posterior
Transverse
LE A Graft 3 Anast
Anterior-Posterior
Transverse
LE A Graft 3 Graft
Anterior-Posterior
Transverse
LE A Graft 3 Native Outflow
Anterior-Posterior
Transverse
LE A Stent 1
Long
Anterior-Posterior
Transverse
LE A Stent 2
Long

Anterior-Posterior
Transverse
LE A Stent 3
Long
Anterior-Posterior
Transverse
LE A Stenosis 1(2D)
Anterior-Posterior
Transverse
Outer Diameter
Inner Diameter
Outer Area
Inner Area
LE A Stenosis 2(2D)
Anterior-Posterior
Transverse
Outer Diameter
Inner Diameter
Outer Area
Inner Area
LE A Stenosis 3(2D)
Anterior-Posterior
Transverse
Outer Diameter
Inner Diameter
Outer Area
Inner Area
LE A Stenosis 4(2D)
Anterior-Posterior
Transverse
Outer Diameter
Inner Diameter
Outer Area
Inner Area
LE A Finding 1
Long
Anterior-Posterior
Transverse
LE A Finding 2
Long
Anterior-Posterior
Transverse
LE A Finding 3
Long
Anterior-Posterior
Transverse
LE A Finding 4
Long

Anterior-Posterior
Transverse
LE A Finding 5
Long
Anterior-Posterior
Transverse
LE A Finding 6
Long
Anterior-Posterior
Transverse
Int Jug V
Anterior-Posterior
Transverse
Checklist
Innom V
Anterior-Posterior
Transverse
Checklist
Subclav V
Anterior-Posterior
Transverse
Checklist
Ax V
Anterior-Posterior
Transverse
Checklist
Brachial V
Anterior-Posterior
Transverse
Checklist
Radial V
Anterior-Posterior
Transverse
Checklist
Ulnar V
Anterior-Posterior
Transverse
Checklist
Volar V
Anterior-Posterior
Transverse
Checklist
Cephalic V
Anterior-Posterior
Transverse
Checklist
Basilic V
Anterior-Posterior

Transverse Checklist	AVF-Outflow Vein Level 1 Anterior-Posterior Transverse
CA Junction Anterior-Posterior Transverse Checklist	AVF-Outflow Vein Level 2 Anterior-Posterior Transverse
Upper Arm Cephalic V Anterior-Posterior Transverse Checklist	AVF-Outflow Vein Level 3 Anterior-Posterior Transverse
Cephalic-Antecubital V Anterior-Posterior Transverse Checklist	AVF-Outflow Vein Level 4 Anterior-Posterior Transverse
Forearm Cephalic V Anterior-Posterior Transverse Checklist	AVF-Outflow Vein Level 5 Anterior-Posterior Transverse
BA Junction Anterior-Posterior Transverse Checklist	AVF-Outflow Vein Level 6 Anterior-Posterior Transverse
Upper Arm Basilic V Anterior-Posterior Transverse Checklist	AVF-Stenosis 1 Anterior-Posterior Transverse
Basilic-Antecubital V Anterior-Posterior Transverse Checklist	AVF-Stenosis 2 Anterior-Posterior Transverse
Forearm Basilic V Anterior-Posterior Transverse Checklist	AVF-Stenosis 3 Anterior-Posterior Transverse
Digital V Anterior-Posterior Transverse Checklist	AVF-Aneurysm 1 Anterior-Posterior Transverse
Median Cubital V Anterior-Posterior Transverse Checklist	AVF-Aneurysm 2 Anterior-Posterior Transverse
AVF-Inflow Artery Anterior-Posterior Transverse	AVF-Aneurysm 3 Anterior-Posterior Transverse
AVF-Anast Anterior-Posterior Transverse	AV Graft-Inflow Artery Anterior-Posterior Transverse
	AV Graft-Arterial Anast Anterior-Posterior Transverse
	AV Graft-Graft Anterior-Posterior Transverse
	AV Graft-Venous Anast Anterior-Posterior Transverse

AV Graft-Outflow Vein Level 1	Transverse
Anterior-Posterior	Checklist
Transverse	Peroneal V
AV Graft-Outflow Vein Level 2	Anterior-Posterior
Anterior-Posterior	Transverse
Transverse	Checklist
AV Graft-Outflow Vein Level 3	Sural V
Anterior-Posterior	Anterior-Posterior
Transverse	Transverse
AV Graft-Outflow Vein Level 4	Checklist
Anterior-Posterior	Soleal V
Transverse	Anterior-Posterior
AV Graft-Outflow Vein Level 5	Transverse
Anterior-Posterior	Checklist
Transverse	A.Tib V
AV Graft-Outflow Vein Level 6	Anterior-Posterior
Anterior-Posterior	Transverse
Transverse	Checklist
C.Iliac V	TP Trunk V
Anterior-Posterior	Anterior-Posterior
Transverse	Transverse
Checklist	Checklist
Ex.Iliac V	Saph V
Anterior-Posterior	Anterior-Posterior
Transverse	Transverse
Checklist	Checklist
IIV	SSV
Anterior-Posterior	Anterior-Posterior
Transverse	Transverse
Checklist	Checklist
CFV	SF Junction
Anterior-Posterior	Anterior-Posterior
Transverse	Transverse
Checklist	Checklist
Femoral V	GSV Thigh
Anterior-Posterior	Anterior-Posterior
Transverse	Transverse
Checklist	Checklist
DFV	GSV Knee
Anterior-Posterior	Anterior-Posterior
Transverse	Transverse
Checklist	Checklist
Pop V	GSV Calf
Anterior-Posterior	Anterior-Posterior
Transverse	Transverse
Checklist	Checklist
P.Tib V	SP Junction
Anterior-Posterior	Anterior-Posterior

Transverse Checklist	P1 PCA
SSV Thigh Extension	P2 PCA
Anterior-Posterior	Ophthalmic A
Transverse	ICA Siphon
Checklist	Terminal Vert A
AASV	BA
Anterior-Posterior	Ba V
Transverse	CCA
Checklist	ICA
PASV	ECA
Anterior-Posterior	Bulb
Transverse	Carotid Bifurcation
Checklist	Vert A
Thigh Perf	Subclav A
Anterior-Posterior	Innom A
Transverse	Mammary A
Checklist	Subclav V
Prox Calf Perf	CCA Aneurysm
Anterior-Posterior	ICA Aneurysm
Transverse	ECA Aneurysm
Checklist	Bulb Aneurysm
Mid Calf Perf	Carotid Bifurcation Aneurysm
Anterior-Posterior	Vert A Aneurysm
Transverse	Subclav A Aneurysm
Checklist	Innom A Aneurysm
Dist Calf Perf	Mammary A Aneurysm
Anterior-Posterior	Carotid Graft 1 Native Inflow
Transverse	Carotid Graft 1 Anast Pre
Checklist	Carotid Graft 1 Anast Max
Pseudoaneurysm	Carotid Graft 1 Anast Post
Long	Carotid Graft 1 Graft
Anterior-Posterior	Carotid Graft 1 Native Outflow
Transverse	Carotid Graft 2 Native Inflow
Neck	Carotid Graft 2 Anast Pre
M-Mode	Carotid Graft 2 Anast Max
D-Mode	Carotid Graft 2 Anast Post
ACA	Carotid Graft 2 Graft
A1 ACA	Carotid Graft 2 Native Outflow
MCA	Carotid Graft 3 Native Inflow
M1 MCA	Carotid Graft 3 Anast Pre
M2 MCA	Carotid Graft 3 Anast Max
ACoM A	Carotid Graft 3 Anast Post
Terminal ICA	Carotid Graft 3 Graft
PCoM A	Carotid Graft 3 Native Outflow
PCA	Carotid Stent 1
	Carotid Stent 2
	Carotid Stent 3
	Axill A

Brachial A
 Ulnar A
 Radial A
 UE A Graft 1 Native Inflow
 UE A Graft 1 Anast
 UE A Graft 1 Graft
 UE A Graft 1 Native Outflow
 UE A Graft 2 Native Inflow
 UE A Graft 2 Anast
 UE A Graft 2 Graft
 UE A Graft 2 Native Outflow
 UE A Graft 3 Native Inflow
 UE A Graft 3 Anast
 UE A Graft 3 Graft
 UE A Graft 3 Native Outflow
 UE A Stent 1
 UE A Stent 2
 UE A Stent 3
 C.Iliac A
 Ex.Iliac A
 IIA
 CFA
 DFA
 SFA
 Pop A
 TP Trunk A
 A.Tib A
 Peroneal A
 P.Tib A
 Dors.Ped. A
 LE A Graft 1 Native Inflow
 LE A Graft 1 Anast Pre
 LE A Graft 1 Anast Max
 LE A Graft 1 Anast Post
 LE A Graft 1 Graft
 LE A Graft 1 Native Outflow
 LE A Graft 2 Native Inflow
 LE A Graft 2 Anast Pre
 LE A Graft 2 Anast Max
 LE A Graft 2 Anast Post
 LE A Graft 2 Graft
 LE A Graft 2 Native Outflow
 LE A Graft 3 Native Inflow
 LE A Graft 3 Anast Pre
 LE A Graft 3 Anast Max
 LE A Graft 3 Anast Post
 LE A Graft 3 Graft
 LE A Graft 3 Native Outflow

LE A Stent 1
 LE A Stent 2
 LE A Stent 3
 Axill V
 Brachial V
 Radial V
 Ulnar V
 Cephalic V
 Basilic V
 AVF-Inflow Artery
 AVF-Anast
 AVF-Outflow Vein Level 1
 AVF-Outflow Vein Level 2
 AVF-Outflow Vein Level 3
 AVF-Outflow Vein Level 4
 AVF-Outflow Vein Level 5
 AVF-Outflow Vein Level 6
 AVF-Stenosis 1
 AVF-Stenosis 2
 AVF-Stenosis 3
 AV Graft-Inflow Artery
 AV Graft-Arterial Anast
 AV Graft-Graft
 AV Graft-Venous Anast
 AV Graft-Outflow Vein Level 1
 AV Graft-Outflow Vein Level 2
 AV Graft-Outflow Vein Level 3
 AV Graft-Outflow Vein Level 4
 AV Graft-Outflow Vein Level 5
 AV Graft-Outflow Vein Level 6
 ASP
 BSP

 CCA(Sten)
 Pre Sten
 Sten
 Post Sten
 ICA(Sten)
 Pre Sten
 Sten
 Post Sten
 ECA(Sten)
 Pre Sten
 Sten
 Post Sten
 Bulb(Sten)
 Pre Sten
 Sten

Post Sten
Carotid Bifurcation(Sten)
Pre Sten
Sten
Post Sten
Vert A(Sten)
Pre Sten
Sten
Post Sten
Subclav A(Sten)
Pre Sten
Sten
Post Sten
Innom A(Sten)
Pre Sten
Sten
Post Sten
Mammary A(Sten)
Pre Sten
Sten
Post Sten
Carotid Stenosis 1
Pre Sten
Sten
Post Sten
Carotid Stenosis 2
Pre Sten
Sten
Post Sten
Carotid Stenosis 3
Pre Sten
Sten
Post Sten
Carotid Stenosis 4
Pre Sten
Sten
Post Sten
Axill A(Sten)
Pre Sten
Sten
Post Sten
Brachial A(Sten)
Pre Sten
Sten
Post Sten
Ulnar A(Sten)
Pre Sten
Sten

Post Sten
Radial A(Sten)
Pre Sten
Sten
Post Sten
UE A Stenosis 1
Pre Sten
Sten
Post Sten
UE A Stenosis 2
Pre Sten
Sten
Post Sten
UE A Stenosis 3
Pre Sten
Sten
Post Sten
UE A Stenosis 4
Pre Sten
Sten
Post Sten
C.Iliac A(Sten)
Pre Sten
Sten
Post Sten
Ex.Iliac A(Sten)
Pre Sten
Sten
Post Sten
IIA(Sten)
Pre Sten
Sten
Post Sten
CFA(Sten)
Pre Sten
Sten
Post Sten
DFA(Sten)
Pre Sten
Sten
Post Sten
SFA(Sten)
Pre Sten
Sten
Post Sten
Pop A(Sten)
Pre Sten
Sten

Post Sten	Checklist
TP Trunk A(Sten)	CFV
Pre Sten	PV
Sten	Reflux
Post Sten	Checklist
A.Tib A(Sten)	Femoral V
Pre Sten	PV
Sten	Reflux
Post Sten	Checklist
Peroneal A(Sten)	DFV
Pre Sten	PV
Sten	Reflux
Post Sten	Checklist
P.Tib A(Sten)	Pop V
Pre Sten	PV
Sten	Reflux
Post Sten	Checklist
Dors.Ped. A(Sten)	P.Tib V
Pre Sten	PV
Sten	Reflux
Post Sten	Checklist
LE A Stenosis 1	Peroneal V
Pre Sten	PV
Sten	Reflux
Post Sten	Checklist
LE A Stenosis 2	Sural V
Pre Sten	PV
Sten	Reflux
Post Sten	Checklist
LE A Stenosis 3	Soleal V
Pre Sten	PV
Sten	Reflux
Post Sten	Checklist
LE A Stenosis 4	A.Tib V
Pre Sten	PV
Sten	Reflux
Post Sten	Checklist
C.Iliac V	TP Trunk V
PV	PV
Reflux	Reflux
Checklist	Checklist
Ex.Iliac V	Saph V
PV	PV
Reflux	Reflux
Checklist	Checklist
IIV	SSV
PV	PV
Reflux	Reflux

Checklist
SF Junction
PV
Reflux
Checklist
GSV Thigh
PV
Reflux
Checklist
GSV Knee
PV
Reflux
Checklist
GSV Calf
PV
Reflux
Checklist
SP Junction
PV
Reflux
Checklist
SSV Thigh Extension
PV
Reflux
Checklist
AASV
PV
Reflux
Checklist
PASV
PV
Reflux
Checklist
Thigh Perf
PV
Reflux
Checklist
Prox Calf Perf
PV
Reflux
Checklist
Mid Calf Perf
PV
Reflux
Checklist
Dist Calf Perf
PV
Reflux

Checklist
ABI
ASP
BSP